

REDGUMbrand
ASSISTED LIVING PRODUCTS

PRODUCT CATALOGUE

Edition 7



E: orders@for-de.com.au (for all email order requests)

P: **08 9248 4180** (Head Office - Perth)

ASK FOR THE BRAND

redgumbrand.com.au

1. **SEAT WALKERS**
2. **SIDE FOLDING WALKERS + ACCESSORIES**
3. **WALKING FRAMES + ACCESSORIES**
4. **WHEELED WALKERS**
5. **KNEE WALKERS**
6. **WALKING STICKS**
7. **WALKING STICKS**
8. **FOLDING WALKING STICKS**
9. **QUAD BASE WALKING STICKS**
10. **WALKING STICK TIPS**
11. **CRUTCHES**
12. **REACHERS**
13. **WHEELCHAIRS**
14. **WHEELCHAIRS + ACCESSORIES**
15. **TRANSIT WHEELCHAIRS**
16. **TILT IN SPACE WHEELCHAIRS**
17. **DRESSING + BATHROOM ACCESSORIES**
18. **SPACE SAVER RANGE**
19. **SPACE SAVER OPTIONS**
20. **BATHROOM PRODUCTS**
21. **MOBILE SHOWER COMMODOES**
22. **BATHROOM + TOILET PRODUCTS**
23. **GRAB RAILS**
24. **KATIE CHAIR RANGE**
25. **KATIE CHAIR ACCESSORIES**
26. **LIFT CHAIRS**
27. **LIFT CHAIRS**
28. **LIFT CHAIRS + ACCESSORIES**
29. **OVER BED CHAIR TABLE RANGE**
30. **AIR CUSHIONS + WEDGE CUSHIONS**
31. **COMFORT FOAM CUSHIONS**
32. **COMFORT LUX BED RANGE**
33. **COMFORT LUX ACCESSORIES**
34. **BED ACCESSORIES**
35. **ALUMINIUM RAMPS**
36. **TELESCOPIC RAMPS**
37. **RUBBER RAMPS**
38. **NOTES**
39. **ABOUT US**



Code: **RG4202AC**
Desc: **ALUM. SEAT WALKER**
PUSH DOWN BRAKE
5" WHEELS
Colours: **RED**
Seat Height: **540mm**
Overall Width: **610mm**
Product Weight: **5.2kg**
Max. Load: **120kg**



Code: **RG4208AJ**
Desc: **ALUM. SEAT WALKER (MINI)**
HAND BRAKE
6" WHEELS
Colours: **RED**
Seat Height: **480mm**
Overall Width: **560mm**
Product Weight: **6kg**
Max. Load: **120kg**



Code: **RG4209**
Desc: **ADJUSTABLE HEIGHT SEAT WALKER**
HAND BRAKE
6" WHEELS
Colours: **PEWTER SILVER / BLUE / RED**
Seat Height: **500 - 560mm**
Overall Width: **580mm**
Product Weight: **6.7kg**
Max. Load: **120kg**



Code: **RG4210**
Desc: **ADJUSTABLE HEIGHT SEAT WALKER**
HAND BRAKE
8" WHEELS
Colours: **PEWTER SILVER / BLUE / RED**
Seat Height: **530 - 590mm**
Overall Width: **600mm**
Product Weight: **7.7kg**
Max. Load: **120kg**



Code: **RG4205**
Desc: **MACK WIDE WALKER**
HAND BRAKE
8" WHEELS
Colours: **RED**
Seat Height: **580mm**
Overall Width: **620mm**
Product Weight: **10kg**
Max. Load: **150kg**



Code: **RG4205A**
Desc: **MINI MACK WIDE WALKER**
HAND BRAKE
8" WHEELS
Colours: **RED**
Seat Height: **500mm**
Overall Width: **625mm**
Product Weight: **9.5kg**
Max. Load: **150kg**



Code: **RG4216C**
Desc: **MIGHTY MACK HD WIDE WALKER**
HAND BRAKE
8" WHEELS
Colours: **RED**
Seat Height: **550mm**
Overall Width: **710mm**
Product Weight: **12.6kg**
Max. Load: **200kg**



RG4212BL / RG4212GY
PRIMO 8 SEAT WALKER
ADJUSTABLE SEAT HEIGHT
8" WHEELS
Colours: **SKY BLUE / GREY**
Seat Height: **505 - 555mm**
Overall Width: **610mm**
Product Weight: **8.3kg**
Max. Load: **136kg**



RG4213GY
PRIMO 8 HEAVY DUTY SEAT WALKER
ADJUSTABLE SEAT HEIGHT
8" WHEELS
Colours: **GREY**
Seat Height: **505 - 555mm**
Overall Width: **690mm**
Product Weight: **9.3kg**
Max. Load: **205kg**





Code: **RG4408**
 Desc: **DUAL PURPOSE SIDE FOLDING WALKER**
 HAND BRAKE, 8" WHEELS
 Colours: **RED / BLUE / WHITE**
 Seat Height: **570mm**
 Overall Width: **600mm**
 Product Weight: **9.3kg**
 Max. Load: **120kg**



Code: **RG4411**
 Desc: **LIGHTWEIGHT SIDE FOLDING WALKER**
 HAND BRAKE, 8" WHEELS
 Colours: **RED / BLUE / WHITE**
 Seat Height: **550mm**
 Overall Width: **615mm**
 Product Weight: **7.5kg**
 Max. Load: **136kg**



Code: **RG4414**
 Desc: **ENDURO SUPER LITE CARBON FIBRE SEAT WALKER**
 HAND BRAKE, 8" WHEELS
 Colours: **BLACK**
 Seat Height: **580mm**
 Overall Width: **610mm**
 Product Weight: **5.9kg**
 Max. Load: **130kg**




Code: **RG4414LS**
 Desc: **ENDURO SUPER LITE CARBON FIBRE SEAT WALKER (LOW SEAT)**
 HAND BRAKE, 8" WHEELS
 Colours: **BLACK**
 Seat Height: **540mm**
 Overall Width: **610mm**
 Product Weight: **5.9kg**
 Max. Load: **130kg**




Code: **RG4401**
 Desc: **COMPACT SIDE FOLDING WALKER**
 8" WHEELS, FOLDS UP
 Colours: **RED / BLUE / WHITE**
 Seat Height: **570mm**
 Overall Width: **640mm**
 Product Weight: **7.5kg**
 Max. Load: **136kg**



Code: **RG4412**
 Desc: **ALL TERRAIN COMPACT SEAT WALKER**
 12" PN TYRES, FOLDS UP
 Colours: **GREY**
 Seat Height: **600mm**
 Overall Width: **610mm**
 Product Weight: **11kg**
 Max. Load: **136kg**




Code: **RG42T**
 Desc: **SEAT WALKER TRAY TABLE**
 (SUITS 500mm WIDE WALKERS)
 Colours: **BLACK**
 Overall Length: **395mm**
 Overall Width: **295mm**
 Product Weight: **0.8kg**
 Max. Load: **8kg**



Code: **RGOBC**
 Desc: **OXYGEN BOTTLE CARRIER**
 SUITABLE FOR WALKERS SUITS 'C' + 'D' SIZE BOTTLE
 Colours: **BLACK**
 Overall Length: **360mm**
 Inner Diameter: **130mm**
 Product Weight: **1.6kg**
 Max. Load: **20kg**



Code: **RG4401CB**
 Desc: **CARRY BAG TO SUIT**
 COMPACT SIDE FOLDING WALKER (RG4401)
 Colours: **DARK BLUE**
 Overall Length: **700mm**
 Overall Width: **400mm**
 Product Weight: **0.5kg**
 Max. Load: **10kg**



Code: **RG040**
 Desc: **WALKING STICK HOLDER**
 TO SUIT MOST FOLDING WALKERS
 Colours: **BLACK**
 Max. Walking Stick Tip Size: **Ø50mm**
 Max. Load: **5kg**



Code: **RG4200**
 Desc: **LOCK AND LOAD CLIP**
 TO KEEP FOLDING SEAT WALKERS CLOSED
 Colours: **RED**



Code: **SE020**
 Desc: **MOBILITY CUP HOLDER**
 ADJUSTABLE CLAMP LOCATE ONTO TUBE LEDGE UP TO Ø40mm
 Colours: **BLACK**
 Bottle Sizes: **55 - 70mm**
 Product Weight: **160g**
 Max. Load: **1kg**



Code: **SE022**
 Desc: **MOBILITY PHONE HOLDER**
 ADJUSTABLE CLAMP LOCATE ONTO TUBE LEDGE UP TO Ø40mm
 Colours: **BLACK**
 Phone Sizes: **60 - 100mm**
 Product Weight: **160g**
 Max. Load: **1kg**



Code: **RG4070MA / RG4070MD**

Desc: **ALUMINIUM WALKER**

ADJUSTABLE HEIGHT FOLDING
FOAM GRIPS

Colours: **SILVER**

Handle Height: **RG4070MA**
750 - 850mm

RG4070MD
815 - 965mm

Internal Width: **500mm**

Product Weight: **4 - 4.5kg**

Max. Load: **150kg**

Code: **RG4070WB**

Desc: **ALUMINIUM WALKER**

ADJUSTABLE HEIGHT FOLDING
5" WHEELS

Colours: **SILVER**

Handle Height: **815 - 965mm**

Internal Width: **500mm**

Product Weight: **4.5kg**

Max. Load: **150kg**

Code: **RG4020AA**

Desc: **PYRAMID ALUM. WALKER**

ADJUSTABLE HEIGHT FOLDING

Colours: **SILVER**

Handle Height: **810 - 910mm**

Internal Width: **500mm**

Product Weight: **3kg**

Max. Load: **150kg**

Code: **RG4020WA**

Desc: **PYRAMID ALUM. WALKER**

ADJUSTABLE HEIGHT FOLDING
5" WHEELS

Colours: **SILVER**

Handle Height: **810 - 910mm**

Internal Width: **500mm**

Product Weight: **3.6kg**

Max. Load: **150kg**



Code: **RG420AG**

Desc: **SWIVEL WHEEL**

Colours: **GREY**

Overall: **345mm (L)**
145mm (W)

Product Weight: **0.5kg**



Code: **RG422AF**

Desc: **WHEEL + GLIDE
BRAKE**

Colours: **GREY**

Overall: **295mm (L)**
100mm (W)

Product Weight: **0.3kg**




Code: **RG410AA**

Desc: **FIXED WHEEL**

Colours: **GREY**

Overall: **240mm (L)**
130mm (W)

Product Weight: **0.3kg**



Code: **RG25**

Desc: **WALKING FRAME
SKI TIP**

TO SUIT Ø25mm
SHAFT
(SOLD IN PAIRS)

Colours: **GREY**

Footprint: **85mm (L)**
40mm (W)

Product Weight: **0.1kg**

Max. Load: **130kg**



Code: **RG25B**

Desc: **WALKING FRAME
SKI TIP**

TO SUIT Ø25mm
SHAFT
(SOLD IN PAIRS)

Colours: **WHITE**

Footprint: **75mm (L)**
30mm (W)

Product Weight: **0.1kg**

Max. Load: **130kg**

® REGISTERED DESIGN



Code: **RG4520K**
 Desc: **TRAY WALKER + SKI TIPS**
 WALKING FRAME & TRANSPORT TROLLEY
 Colours: **SILVER**
 Overall Height: **820 - 945mm**
 Internal Width: **495mm**
 Product Weight: **8kg**
 Max. Load: **150kg**

Code: **RG4300AA**
 Desc: **WALKING TUTOR**
 ADJUSTABLE HEIGHT
 8" WHEELS, SEAT
 Colours: **GREY**
 Armrest Height: **1008 - 1296mm**
 Overall Width: **900mm**
 Product Weight: **11kg**
 Max. Load: **140kg**

Code: **RG4310AA**
 Desc: **WALKING TUTOR**
 ADJUSTABLE HEIGHT
 8" WHEELS, SEAT
 Colours: **RED**
 Armrest Height: **975 - 1115mm**
 Overall Width: **610mm**
 Product Weight: **14.2kg**
 Max. Load: **120kg**

Code: **RG4402**
 Desc: **TRI WALKER**
 ADJUSTABLE HEIGHT FOLDING
 8" WHEELS
 Colours: **CHAMPAGNE SKY BLUE**
 Hand Height: **830 - 930mm**
 Overall Width: **675mm**
 Product Weight: **4.6kg**
 Max. Load: **100kg**

Code: **RG4404**
 Desc: **TRI WALKER**
 ADJUSTABLE HEIGHT FOLDING
 8" WHEELS
 Colours: **RED / BLUE**
 Hand Height: **875 - 1025mm**
 Overall Width: **610mm**
 Product Weight: **7.5kg**
 Max. Load: **100kg**



Code: **RG20KWB**
Desc: **KNEE WALKER**
ADJUSTABLE HEIGHT
HAND BRAKE, 8" WHEELS
Colours: **BLUE**
Knee Platform Height: **410 - 510mm**
Handlebar Height: **870 - 1080mm**
Product Weight: **38kg**
Max. Load: **135kg**

Code: **RG25KWB**
Desc: **COMFORT-RIDE**
KNEE WALKER
12 / 10" PNEUMATIC TYRES
Colours: **METALLIC GREY**
Knee Platform Height: **470 - 600mm**
Handlebar Height: **850 - 1050mm**
Product Weight: **13kg**
Max. Load: **136kg**

Code: **RG30KWB**
Desc: **KNEE WALKER**
ADJUSTABLE HEIGHT
HAND BRAKE, 8" WHEELS
Colours: **COPPER**
Knee Platform Height: **425 - 550mm**
Handlebar Height: **800 - 940mm**
Product Weight: **12.1kg**
Max. Load: **130kg**

Code: **RG40KWH**
Desc: **HEAVY DUTY**
KNEE WALKER
ADJUSTABLE HEIGHT
HAND BRAKE, 8" WHEELS
Colours: **COPPER**
Knee Platform Height: **425 - 550mm**
Handlebar Height: **800 - 940mm**
Product Weight: **14.5kg**
Max. Load: **180kg**





Code: **RGWS320 / RGWS321
RGWS322**
Desc: **HYBRID WALKING STICK
POLE**
ADJUSTABLE HEIGHT
Colours: **BLACK CAMO
BLUE CAMO
GREEN CAMO**
Overall Height: **930- 1200mm**
Product Weight: **0.5kg**
Max. Load: **113kg**



Code: **RGWS330 / RGWS331**
Desc: **RISE WALKING STICK
ADJUSTABLE HEIGHT**
Colours: **SMOKE GREY / BLUE**
Overall Height: **700mm**
Product Weight: **0.4kg**
Max. Load: **113kg**



Code: **RGWS340 / RGWS341
RGWS342**
Desc: **SOFT TOUCH NARROW
NECK WALKING STICK**
METAL JOINT CONNECTION
Colours: **GEO BLACK
GEO ORANGE
GEO DARK BLUE**
Overall Height: **760 - 990mm**
Product Weight: **0.4kg**
Max. Load: **113kg**




Code: **RGWS300**
Desc: **FOLDING HIKING
WALKING STICK**
PERFECT FOR HIKING
Colours: **BROWN**
Overall Height: **1260 -1400mm**
Product Weight: **0.5kg**
Max. Load: **113kg**



Code: **RG1040**
Desc: **WALKING STICK
OFFSET CURVED
HANDLE**
ADJUSTABLE HEIGHT
Colours: **BLACK / BROWN**
Overall Height: **765 - 990mm**
Product Weight: **0.35kg**
Max. Load: **100kg**



Code: **RG1060**
Desc: **WALKING STICK
OFFSET STRAIGHT
HANDLE**
ADJUSTABLE HEIGHT
Colours: **BLACK / BROWN**
Overall Height: **765 - 990mm**
Product Weight: **0.35kg**
Max. Load: **100kg**



Code: **RGWS351B
RGWS351G**
Desc: **WALKING STICK L/H
ANATOMIC HANDLE**
ADJUSTABLE HEIGHT
Colours: **BLUE / GREY**
Overall Height: **760 - 990mm**
Product Weight: **0.4kg**
Max. Load: **113kg**



Code: **RGWS350B
RGWS350G**
Desc: **WALKING STICK R/H
ANATOMIC HANDLE**
ADJUSTABLE HEIGHT
Colours: **BLUE / GREY**
Overall Height: **760 - 990mm**
Product Weight: **0.4kg**
Max. Load: **113kg**



Code: **RG2060**
Desc: **WALKING STICK "T"
PLASTIC HANDLE**
ADJUSTABLE HEIGHT
Colours: **BLACK**
Overall Height: **765 - 990mm**
Product Weight: **0.3kg**
Max. Load: **120kg**



Code: **RG2010**
Desc: **WALKING STICK "T"
WOODEN HANDLE**
ADJUSTABLE HEIGHT
Colours: **BLACK / BROWN**
Overall Height: **765 - 990mm**
Product Weight: **0.3kg**
Max. Load: **120kg**



Code: **RG2010A**
Desc: **WALKING STICK
SWAN WOODEN NECK**
ADJUSTABLE HEIGHT
Colours: **BLACK / BROWN**
Overall Height: **765 - 990mm**
Product Weight: **0.35kg**
Max. Load: **120kg**



Code: **RG2015**
Desc: **WALKING STICK
COMFORT GRIP**
ADJUSTABLE HEIGHT
Ø19mm SHAFT
Colours: **RED**
Overall Height: **710 - 950mm**
Product Weight: **0.35kg**
Max. Load: **120kg**



Code: **RG2025**
Desc: **WALKING STICK
COMFORT GRIP**
ADJUSTABLE HEIGHT
Ø15mm SHAFT
Colours: **BLACK**
Overall Height: **710 - 950mm**
Product Weight: **0.35kg**
Max. Load: **120kg**



Code: **RGWS201**
Desc: **CROOKED CANE
WITH HAND STRAP**
ADJUSTABLE HEIGHT,
ASSIST USER GETTING
UP FROM SEATED
POSITION
Colours: **BRONZE**
Overall Height: **810 - 1040mm**
Product Weight: **0.75kg**
Max. Load: **113kg**



Code: **RGWSH**
Desc: **WALKING STICK
HOLDER**
CLIPS ONTO ANY
WALKING STICK
WITH 22mm OR LESS
Colours: **BLACK/ ORANGE**
Product Weight: **0.3kg**




WALKING STICK HOLDER



- ORANGE REFLECTIVE STRIP THAT ACTS AS A VISUAL AID WHEN WALKING
- CLIPS ONTO MOST WALKING STICKS
- EASY TO FIT AND REMOVE
- SPRING LOADED LEDGE WITH FELT PAD ON UNDERSIDE
- ALLOWS THE WALKING STICK TO BE HELD UPRIGHT ON A TABLE CUPBOARD



<p>Code: RGWS310 Desc: WALKING STICK PREMIUM TELESCOPIC ADJUSTABLE HEIGHT ALUMINIUM Colours: PAISLEY Overall Height: 750 - 910mm Product Weight: 0.3kg Max. Load: 113kg</p>	<p>Code: RGWS311 Desc: WALKING STICK PREMIUM TELESCOPIC ADJUSTABLE HEIGHT ALUMINIUM Colours: DIAMOND Overall Height: 750 - 910mm Product Weight: 0.3kg Max. Load: 113kg</p>	<p>Code: RGWS312 Desc: WALKING STICK PREMIUM TELESCOPIC ADJUSTABLE HEIGHT ALUMINIUM Colours: FLOWERS Overall Height: 750 - 910mm Product Weight: 0.3kg Max. Load: 113kg</p>	<p>Code: RG3020 Desc: WALKING STICK FOLDING "T" HANDLE ADJUSTABLE HEIGHT ALUMINIUM Colours: BLACK SILVER Overall Height: 840 - 940mm Product Weight: 0.4kg Max. Load: 120kg</p>	<p>Code: RG3020/26 Desc: WALKING STICK FOLDING "T" HANDLE ADJUSTABLE HEIGHT ALUMINIUM Colours: BLACK ROSE Overall Height: 840 - 940mm Product Weight: 0.4kg Max. Load: 120kg</p>	<p>Code: RG3020/01 Desc: WALKING STICK FOLDING "T" HANDLE ADJUSTABLE HEIGHT ALUMINIUM Colours: BLUE CHECK Overall Height: 840 - 940mm Product Weight: 0.4kg Max. Load: 120kg</p>	<p>Code: RG3020/19 Desc: WALKING STICK FOLDING "T" HANDLE ADJUSTABLE HEIGHT ALUMINIUM Colours: SUMMER FLORAL Overall Height: 840 - 940mm Product Weight: 0.4kg Max. Load: 120kg</p>	<p>Code: RG3020/20 Desc: WALKING STICK FOLDING "T" HANDLE ADJUSTABLE HEIGHT ALUMINIUM Colours: BLUE FLORAL Overall Height: 840 - 940mm Product Weight: 0.4kg Max. Load: 120kg</p>	<p>Code: RG3020/29 Desc: WALKING STICK FOLDING "T" HANDLE ADJUSTABLE HEIGHT ALUMINIUM Colours: ARCTIC FLORAL Overall Height: 840 - 940mm Product Weight: 0.4kg Max. Load: 120kg</p>
---	---	---	---	--	--	---	---	---



<p>Code: RGWS301 / RGWS302 / RGWS303 / RGWS304 / RGWS305 / RGWS306 Desc: FOLDING CANE ADJUSTABLE HEIGHT ALUMINIUM HAND STRAP Colours: RED "PEARL" / OLIVE "PEARL" / DARK BLUE "PEARL" / SILVER GREY "PEARL" / AQUA BLUE "PEARL" / BLACK "PEARL" Overall Height: 825 - 925mm Product Weight: 0.4kg Max. Load: 113kg</p>	<p>Code: RGWS307 Desc: PREMIUM FOLDING WALKING STICK ALUMINIUM FRAME ADJUSTABLE HEIGHT Colours: SPARKLE BLACK Overall Height: 830 - 930mm Product Weight: 0.35kg Max. Load: 113kg</p>	<p>Code: RGWS308 Desc: PREMIUM FOLDING WALKING STICK ALUMINIUM FRAME ADJUSTABLE HEIGHT Colours: SPARKLE GOLD Overall Height: 830 - 930mm Product Weight: 0.35kg Max. Load: 113kg</p>	<p>Code: RGWS309 Desc: PREMIUM FOLDING WALKING STICK ALUMINIUM FRAME ADJUSTABLE HEIGHT Colours: SPARKLE BLUE Overall Height: 830 - 930mm Product Weight: 0.35kg Max. Load: 113kg</p>	<p>Code: RG3091 Desc: TRIPOD CANE FOLDABLE WITH SEAT ALUMINIUM FRAME PLASTIC SEAT Colours: GREEN Seat Height: 545mm Product Weight: 0.8kg Max. Load: 100kg</p>
--	---	--	--	--



Code: **RG5610**
 Desc: **STEEL QUAD CANE
 OFFSET HANDLE**
 ADJUSTABLE HEIGHT
 310 x 205 x 75mm
 BASE
 Colours: **SILVER / BLACK**
 Overall Height: **755 - 980mm**
 Product Weight: **1kg**
 Max. Load: **100kg**

Code: **RG5620**
 Desc: **STEEL QUAD CANE
 OFFSET HANDLE**
 ADJUSTABLE HEIGHT
 230 x 155 x 75mm
 BASE
 Colours: **SILVER / BLACK**
 Overall Height: **755 - 980mm**
 Product Weight: **0.9kg**
 Max. Load: **100kg**

Code: **RG5650**
 Desc: **ALUM. QUAD CANE
 OFFSET HANDLE**
 ADJUSTABLE HEIGHT
 305 x 240 x 250mm
 BASE
 Colours: **COPPER**
 Overall Height: **790 - 995mm**
 Product Weight: **0.7kg**
 Max. Load: **100kg**

Code: **RGWS101**
 Desc: **CROOKED QUAD
 CANE WITH
 HAND STRAP**
 ADJUSTABLE HEIGHT
 ASSIST USER
 GETTING UP FROM
 SEATED POSITION
 Colours: **BRONZE**
 Overall Height: **850 - 1080mm**
 Product Weight: **0.9kg**
 Max. Load: **113kg**

Code: **RG5611**
 Desc: **HEAVY DUTY
 QUAD CANE WITH
 OFFSET HANDLE**
 ADJUSTABLE HEIGHT
 300 x 210 x 80mm
 BASE
 Colours: **BLACK**
 Overall Height: **740 - 970mm**
 Product Weight: **1.7kg**
 Max. Load: **227kg**






Code: **RGT103**
 Desc: **SELF STANDING WALKING STICK TIP**
 TO SUIT Ø19mm SHAFT
 Colours: **GREY**
 Footprint: **Ø115mm**
 Max. Load: **130kg**

Code: **RGT102 BLACK**
 Desc: **SELF STANDING WALKING STICK TIP**
 TO SUIT Ø19mm SHAFT
 Colours: **BLACK**
 Footprint: **Ø115mm**
 Max. Load: **130kg**

Code: **RGT101**
 Desc: **SELF STANDING WALKING STICK TIP**
 FLEXIBLE BASE TO SUIT Ø16 / 19mm SHAFT
 Colours: **BLACK**
 Footprint: **Ø115mm**
 Max. Load: **130kg**



Code: **RG16**
 Desc: **WALKING STICK TIP**
 TO SUIT Ø16mm SHAFT
 Colours: **BLACK / GREY**
 Footprint: **Ø30mm**
 Max. Load: **130kg**

Code: **RG19**
 Desc: **WALKING STICK TIP**
 TO SUIT Ø19mm SHAFT
 Colours: **BLACK / GREY**
 Footprint: **Ø37.5mm**
 Max. Load: **130kg**

Code: **RG21**
 Desc: **WALKING STICK TIP**
 TO SUIT Ø22mm SHAFT
 Colours: **WHITE**
 Footprint: **Ø50mm**
 Max. Load: **130kg**

Code: **RG22**
 Desc: **WALKING STICK TIP**
 TO SUIT Ø22mm SHAFT
 Colours: **GREY**
 Footprint: **Ø45mm**
 Max. Load: **130kg**

Code: **RG23**
 Desc: **WALKING STICK TIP**
 TO SUIT Ø28.6mm SHAFT
 Colours: **GREY**
 Footprint: **Ø40mm**
 Max. Load: **130kg**

Code: **RG26**
 Desc: **WALKING STICK TIP**
 TO SUIT Ø25.4mm SHAFT
 Colours: **GREY**
 Footprint: **Ø57.5mm**
 Max. Load: **130kg**

Code: **RG27**
 Desc: **WALKING STICK TIP**
 TO SUIT Ø25.4mm SHAFT
 Colours: **GREY**
 Footprint: **Ø40mm**
 Max. Load: **130kg**

Code: **RG28**
 Desc: **WALKING STICK TIP**
 TO SUIT Ø28.6mm SHAFT
 Colours: **GREY**
 Footprint: **Ø40mm**
 Max. Load: **130kg**



Code: **RG7801**
 Desc: **ALUMINIUM FOREARM CRUTCH**
 ANATOMIC GRIPS (L/R)
 HEIGHT ADJUSTABLE
 PLASTIC SILENCER
 Colours: **GREY**
 Overall Height: **600 - 830mm**
 Product Weight: **0.8kg (each)**
 Max. Load: **120kg (pair)**



Code: **RG7802**
 Desc: **ALUMINIUM FOREARM CRUTCH**
 ANATOMIC GRIPS (L/R)
 HEIGHT ADJUSTABLE
 PLASTIC SILENCER
 Colours: **GREY**
 Overall Height: **670 - 900mm**
 Product Weight: **0.8kg (each)**
 Max. Load: **120kg (pair)**

Code: **RG7803**
 Desc: **ALUMINIUM FOREARM CRUTCH**
 ANATOMIC GRIPS (L/R)
 HEIGHT ADJUSTABLE
 PLASTIC SILENCER
 Colours: **GREY**
 Overall Height: **780 - 1060mm**
 Product Weight: **0.8kg (each)**
 Max. Load: **120kg (pair)**






Code: **RG7804**
 Desc: **ALUMINIUM FOREARM CRUTCH**
 UNIVERSAL GRIP HEIGHT
 ADJUSTABLE REFLECTOR
 Colours: **GREY**
 Overall Height: **780 - 1060mm**
 Product Weight: **0.6kg (each)**
 Max. Load: **113kg (pair)**

NEW!




Code: **RG7200M / RG7200T**
 Desc: **ADJUSTABLE UNDERARM CRUTCHES**
 UNIVERSAL GRIP HEIGHT
 ADJUSTABLE REFLECTOR
 Colours: **GREY**
 Overall Height: **1150 - 1350mm (RG7200M)**
 1350 - 1550mm (RG7200T)
 Product Weight: **1.6kg (each) (RG7200M)**
 1.8kg (each) (RG7200T)
 Max. Load: **136kg**



Code: **RGM24**
 Desc: **DELUXE REACHER 24"**
 ALUMINIUM BODY, MAGNET + HOOK
 Colours: **SILVER WITH BLACK HANDLE**
 Overall Length: **600mm**
 Product Weight: **0.17kg**
 Max. Load: **1.5kg**

Code: **RGM30**
 Desc: **DELUXE REACHER 30"**
 ALUMINIUM BODY, MAGNET + HOOK
 Colours: **SILVER WITH BLACK HANDLE**
 Overall Length: **760mm**
 Product Weight: **0.18kg**
 Max. Load: **1.5kg**

Code: **RGM35**
 Desc: **DELUXE REACHER 35"**
 ALUMINIUM BODY, MAGNET + HOOK
 Colours: **SILVER WITH BLACK HANDLE**
 Overall Length: **890mm**
 Product Weight: **0.2kg**
 Max. Load: **1.5kg**



Code: **RGM32**
 Desc: **REACHER 32"**
 ALUMINIUM BODY, 90° ROTATABLE HEAD
 Colours: **SILVER WITH BLACK HANDLE**
 Overall Length: **820mm**
 Product Weight: **0.2kg**
 Max. Load: **1.5kg**

Code: **RGMF32**
 Desc: **FOLDING REACHER 32"**
 ALUMINIUM BODY, CUPPED CLAW
 LOCKING MECHANISM
 Colours: **SILVER WITH BLUE HANDLE**
 Overall Length: **820mm**
 Product Weight: **0.2kg**
 Max. Load: **1.5kg**

Code: **RGM36HD**
 Desc: **HEAVY DUTY REACHER**
 ALUMINIUM BODY, ROTATABLE HEAD
 WITH LOCKING MECHANISM
 Colours: **SILVER WITH BLACK HANDLE**
 Overall Length: **870mm**
 Product Weight: **0.4kg**
 Max. Load: **5kg**



Code: **RG310A18BK / RG310A18BL / RG310A18R**

Desc: **QUARTZ WHEELCHAIR**

ALUMINIUM, LIFT UP ARMS, PU TYRES, SWING AWAY LEG RESTS, QUICK RELEASE RIMS, FOLDS FOR STORAGE

Colours: **BLACK / BLUE / RED**

Seat Width: **19" (485mm)**

Product Weight: **13.5kg**

Max. Load: **115kg**



Code: **RG314GR18 / RG314BL18**

Desc: **COMFORTLITE WHEELCHAIR**

ALUMINIUM, PU TYRES, SWING AWAY LEG RESTS, FOLDS FOR STORAGE

Colour: **GREEN / STEEL BLUE**

Seat Width: **18" / 460mm**

Product Weight: **12kg**

Max. Load: **100kg**



Code: **RG18B**

Desc: **CRYSTAL WHEELCHAIR**

ALUMINIUM, SOLID PU TYRES, SWING AWAY LEG RESTS, QUICK RELEASE RIMS FOLDS FOR STORAGE

Colours: **SILVER**

Seat Width: **450mm**

Product Weight: **12.5kg**

Max. Load: **127kg**



Code: **RG20B**

Desc: **CRYSTAL WHEELCHAIR**

ALUMINIUM, SOLID PU TYRES, SWING AWAY LEG RESTS, QUICK RELEASE RIMS FOLDS FOR STORAGE

Colours: **SILVER**

Seat Width: **490mm**

Product Weight: **12.5kg**

Max. Load: **127kg**



FOLDING FRAME

FOLD DOWN BACKREST

RG310A18BK / RG310A18BL / RG310A18R
QUARTZ WHEELCHAIR



FOLDS FOR STORAGE AND TRANSIT

RG314GR18 / RG314BL18
COMFORTLITE WHEELCHAIR



RG18B / RG20B
CRYSTAL WHEELCHAIR



Code: **RG18AAD**
Desc: **PEARL 18" WHEELCHAIR**
ALUMINIUM, LIFT UP ARMS, SOLID PU TYRES, SWING AWAY LEG RESTS, QUICK RELEASE RIMS
Colours: **PEARL WHITE**
Seat Width: **18" (460mm)**
Product Weight: **21.1kg**
Max. Load: **120kg**



Code: **RG20AAD**
Desc: **PEARL 20" WHEELCHAIR**
ALUMINIUM, LIFT UP ARMS, SOLID PU TYRES, SWING AWAY LEG RESTS, QUICK RELEASE RIMS
Colours: **PEARL WHITE**
Seat Width: **20" (510mm)**
Product Weight: **22.5kg**
Max. Load: **130kg**



Code: **RG18AADRB**
Desc: **PEARL 18" WHEELCHAIR REAR HAND BRAKE**
OTHERWISE SAME SPECS AS RG18AAD: PEARL WHEELCHAIR
Colours: **PEARL WHITE**
Seat Width: **18" (460mm)**
Product Weight: **22kg**
Max. Load: **120kg**



Code: **RG20AADRB**
Desc: **PEARL 20" WHEELCHAIR REAR HAND BRAKE**
OTHERWISE SAME SPECS AS RG20AAD: PEARL WHEELCHAIR
Colours: **PEARL WHITE**
Seat Width: **20" (510mm)**
Product Weight: **23kg**
Max. Load: **130kg**



Code: **RG22HD**
Desc: **TITAN WHEELCHAIR**
SUPER STRONG STEEL FRAME, ADJUSTABLE BACKREST, LIFT UP ARMS, SWING LEG RESTS
Colours: **MARBLE BLACK**
Seat Width: **22" (560mm)**
Product Weight: **25kg**
Max. Load: **200kg**



Code: **RG18TW**
Desc: **TRAKKA WHEELCHAIR**
ULTRA STRONG STEEL FRAME 45mm WIDE PU TYRES SWING AWAY LEG RESTS
Colours: **GUNMETAL GREY**
Seat Width: **18" (460mm)**
Product Weight: **19.5kg**
Max. Load: **150kg**



Code: **RG20TW**
Desc: **TRAKKA WHEELCHAIR**
ULTRA STRONG STEEL FRAME 45mm WIDE PU TYRES SWING AWAY LEG RESTS
Colours: **GUNMETAL GREY**
Seat Width: **20" (510mm)**
Product Weight: **21.5kg**
Max. Load: **150kg**



Code: **RGCWCC**
Desc: **COMPACT WHEELCHAIR CARRIER**
FITS TO MOST TOW BARS
Colours: **GOLD ZINC PLATED**
Internal Length: **535mm**
Internal Depth: **350mm**
Product Weight: **10.5kg**
Max. Load: **32kg**

® REGISTERED DESIGN



Code: **RGCWCCB**
Desc: **COMPACT WHEELCHAIR CARRIER CARRIER BASKET INSERT**
TO SUIT TRANSIT WHEELCHAIRS
Colours: **GOLD ZINC PLATED**
Internal Length: **535mm**
Internal Depth: **350mm**
Product Weight: **3kg**
Max. Load: **32kg**



Code: **RGIVOCB**
Desc: **I.V POLE + OXYGEN BOTTLE CARRY FRAME**
SUITS POWER CHAIR UNIT ONLY
Colours: **CHROME**
Inner Diameter: **Ø130mm**
Product Weight: **2kg**
Max. Load: **20kg**

NEW!



Code: **650069**
Desc: **WHEELCHAIR REPLACEMENT TYRE**
REPLACEMENT TYRE FOR SELF PROPELLED WHEELCHAIRS
Colours: **GREY**
Overall: **Ø24" x 1-3/8"**
Product Weight: **0.6kg**

ADDED COMFORT



Refer to pages 30 -31 for information on CUSHION RANGE



Code: **RG312GR18 / RG312L18**
 Desc: **COMFORTLITE TRANSIT WHEELCHAIR**
 ALUMINIUM, FOLDING BACKREST, PU TYRES, ANTI TIPPERS
 Colour: **GREEN / STEEL BLUE**
 Seat Width: **18" (460mm)**
 Product Weight: **9.5kg**
 Max. Load: **100kg**



Code: **RG311Y18 / RG311B18**
 Desc: **OPAL TRANSIT WHEELCHAIR**
 ALUMINIUM, FOLDING BACKREST, PU TYRES, ANTI TIPPERS
 Colour: **YELLOW / BLUE**
 Seat Width: **18" (460mm)**
 Product Weight: **9.4kg**
 Max. Load: **136kg**



Code: **RG17HBD**
 Desc: **17" TRANSPORT CHAIR DELUXE**
 FOLDING BACKREST, HAND BRAKES, REMOVABLE ARMRESTS, SWING AWAY LEG RESTS
 Colours: **RED/BLUE**
 Seat Width: **17" (430mm)**
 Product Weight: **10.5kg**
 Max. Load: **130kg**



Code: **RG19HBD**
 Desc: **19" TRANSPORT CHAIR DELUXE**
 FOLDING BACKREST, HAND BRAKES, REMOVABLE ARMRESTS, SWING AWAY LEG RESTS
 Colours: **RED / BLUE**
 Seat Width: **19" (485mm)**
 Product Weight: **10.5kg**
 Max. Load: **130kg**



Code: **RG22HBD**
 Desc: **22" TRANSPORT CHAIR DELUXE**
 FOLDING BACKREST, HAND BRAKES, REMOVABLE ARMRESTS, SWING AWAY LEG RESTS
 Colours: **RED**
 Seat Width: **22" (560mm)**
 Product Weight: **13kg**
 Max. Load: **180kg**



FOLDS UP FOR TRANSIT AND STORAGE





Code: **RG350TR18 / RG350TR20**

Desc: **ONYX TILT IN SPACE
TRANSIT WHEELCHAIR
450mm / 500mm**

HEIGHT ADJUSTABLE SEAT
AND BACKREST

Colours: **BLACK**

Overall Height: **990 - 1060mm**

Overall Width: **680mm (450mm)
730mm (500mm)**

Seat Width: **350 - 500mm (450mm)
400 - 550mm (500mm)**

Seat Depth: **460mm - 520mm**

Product Weight: **38kg**

Max. Load: **135kg**



Code: **RG350SP18 / RG350SP20**

Desc: **ONYX TILT IN SPACE
SELF PROPELLED
WHEELCHAIR
450mm / 500mm**

HEIGHT ADJUSTABLE SEAT
AND BACKREST

Colours: **BLACK**

Overall Height: **990 - 1060mm**

Overall Width: **680mm (450mm)
730mm (500mm)**

Seat Width: **350 - 500mm (450mm)
400 - 550mm (500mm)**

Seat Depth: **460 - 520mm**

Product Weight: **38kg**

Max. Load: **135kg**

These wheelchairs have been crash tested according to standard "ISO 7176-19" and may be used as a seat in motor vehicles.




**AVAILABLE IN
18" AND 20"
SEAT WIDTH**




**QUICK RELEASE
REAR WHEELS**

® REGISTERED DESIGN




NEW!



NEW!



NEW!



NEW!



Code: **RG71SH**
Desc: **SHOE HORN DRESSING STICK**
IDEAL FOR PEOPLE WITH LIMITED REACH
Colours: **RED / BLUE / SILVER**
Overall Length: **655mm**
Product Weight: **0.1kg**


Code: **RG71DS**
Desc: **DRESSING STICK**
ASSIST WITH DRESSING AND REMOVING CLOTHING
Colours: **WOOD**
Overall: **690mm (L)**
Handle: **Ø15mm**
Hook: **Ø6 / Ø3mm**
Product Weight: **0.07kg**

Code: **RG71LL**
Desc: **LEG LIFT**
ASSIST WITH RAISING AND LOWERING LEG
Colours: **BLUE**
Overall: **950mm (L)**
Strap: **25mm (W)**
Loop Handle: **160mm (L)**
Product Weight: **0.01kg**


Code: **RG71SA**
Desc: **SOCK AID**
ASSIST WITH PUTTING SOCKS ON WITH NO BENDING
Colours: **BLUE / WHITE**
Overall: **875mm (L) x 180mm (W)**
3mm (thick)
Sock Holder: **225mm (L)**
Strap: **650 (L) X 25mm (W)**
Loop Handle: **105mm (L)**
Product Weight: **0.09kg**

Code: **RG71MUB**
Desc: **MALE URINAL BOTTLE**
ASSIST MEN WITH TOILETING AWAY FROM TOILET
Colours: **WHITE**
Overall: **225mm (L) x 95mm (W)**
130mm (H)
Sprout int: **Ø52mm**
Capacity: **1 litre**
Product Weight: **0.12kg**


NEW!




NEW!




NEW!



NEW!



NEW!



Code: **RG71HB**
Desc: **LONG HANDLE HAIR BRUSH**
ASSIST WITH LONGER REACH FOR HAIR BRUSHING
Colours: **BLUE/WHITE**
Overall: **355mm (L)**
Brush: **105mm (L) x 50mm (W)**
25mm (H)
Product Weight: **0.11kg**


Code: **RG71HC**
Desc: **LONG HANDLE HAIR COMB**
ASSIST WITH LONGER REACH FOR HAIR COMBING
Colours: **BLUE/WHITE**
Overall: **380mm (L)**
Comb: **100mm (L) x 50mm (W)**
27mm (H)
Product Weight: **0.08kg**

Code: **RG71HS**
Desc: **LONG HANDLE SPONGE**
IDEAL FOR WASHING THOSE HARD TO REACH PLACES
Colours: **BLUE/WHITE**
Overall: **650mm (L)**
Sponge: **125mm (L) x 75mm (W)**
50mm (H)
Handle: **525mm (L) x Ø12mm**
Product Weight: **0.1kg**

Code: **RG71RBM**
Desc: **RECTANGULAR BATH MAT**
PROVIDES A SECURE NON-SLIP SURFACE IN A WET AREA
Colours: **GREY / BLUE**
Size: **765mm (L) x 345mm (W)**
5mm (H)
Product Weight: **0.68kg**

Code: **RG71SBM**
Desc: **SQUARE BATH MAT**
PROVIDES A SECURE NON-SLIP SURFACE IN A WET AREA
Colours: **GREY / BLUE**
Size: **550mm (L) x 550mm (W)**
5mm (H)
Product Weight: **0.68kg**

® REGISTERED DESIGN



® PATENTED DESIGN

IMPROVED DESIGN



® PATENTED DESIGN



® PATENTED DESIGN



COMPLIANT TO RELEVANT AUST / NZ STANDARDS

COMPLIANT TO RELEVANT AUST / NZ STANDARDS

COMPLIANT TO RELEVANT AUST / NZ STANDARDS


Code: **RG553**
 Desc: **SPACE SAVER CORNER SHOWER STOOL**
 ADJUSTABLE HEIGHT,
 FLAT PACK CONTOURED SEAT
 Colours: **WHITE with SILVER LEGS**
 Seat Height: **490 - 615mm**
 Seat Size: **395mm (W)**
295mm (D)
 Product Weight: **1.5kg**
 Max. Load: **130kg**

Code: **RG554H2**
 Desc: **SPACE SAVER SHOWER STOOL**
 ADJUSTABLE HEIGHT
 FLAT PACK
 Colours: **WHITE with SILVER LEGS**
 Seat Height: **485 - 610mm**
 Seat Size: **470mm (W)**
355mm (D)
 Product Weight: **2.6kg**
 Max. Load: **130kg**


Code: **RG555H**
 Desc: **SPACE SAVER SHOWER CHAIR**
 ADJUSTABLE HEIGHT,
 FLAT PACK CONTOURED BACKREST
 Colours: **WHITE with SILVER LEGS**
 Seat Height: **485 - 610mm**
 Seat Size: **470mm (W)**
355mm (D)
 Product Weight: **3.1kg**
 Max. Load: **130kg**

Code: **RG515**
 Desc: **SPACE SAVER TOILET SEAT RAISER**
 ADJUSTABLE HEIGHT,
 FLAT PACK, SPLASHGUARD
 Colours: **WHITE with SILVER LEGS**
 Seat Height: **500 - 625mm**
 Seat Size: **470mm (W)**
400mm (D)
 Product Weight: **3.5kg**
 Max. Load: **130kg**


DEMENTIA ASSIST



Code: **RG5SHR / RG5SHB**
 Desc: **SPACE SAVER DEMENTIA ASSIST HANDLES**
 DESIGNED TO ASSIST PEOPLE SUFFERING WITH DEMENTIA
 Colour: **RED / BLUE**
 Product Weight: **0.3kg**
 Max. Load: **130kg**



Code: **RG53GP**
 Desc: **SPACE SAVER CORNER SHOWER STOOL GEL PAD**
 10mm GEL PAD TO RELIEVE PRESSURE WHEN SITTING
 Colour: **BLUE**
 Overall: **365mm (W) x 265mm (D) 10mm (H)**
 Product Weight: **1kg**



Code: **RG54GP**
 Desc: **SPACE SAVER SHOWER STOOL GEL PAD**
 10mm GEL PAD TO RELIEVE PRESSURE WHEN SITTING
 Colour: **BLUE**
 Overall: **355mm (W) x 285mm (D) 10mm (H)**
 Product Weight: **1kg**



Code: **RG55GP**
 Desc: **SPACE SAVER SHOWER CHAIR GEL PAD**
 10mm GEL PAD TO RELIEVE PRESSURE WHEN SITTING
 Colour: **BLUE**
 Overall: **312mm (W) x 350mm (D) 10mm (H)**
 Product Weight: **1kg**




Code: **RG51GP**
 Desc: **SPACE SAVER TOILET SEAT RAISER GEL PAD**
 10mm GEL PAD TO RELIEVE PRESSURE WHEN SITTING
 Colour: **BLUE**
 Overall: **370mm (W) x 380mm (D) 10mm (H)**
 Product Weight: **1kg**


REGISTERED DESIGN



Code: **RG55LH**
 Desc: **SPACE SAVER LOW HANDLES**
 MORE SPACE TO OPERATE WITH EASE
 Colour: **WHITE**
 Overall: **190mm (W) x 20mm (D) 110mm (H)**
 Product Weight: **0.3kg**
 Max. Load: **130kg**



Code: **RG55C**
 Desc: **SPACE SAVER SHOWER CADDY**
 IDEAL FOR ASSISTING AND ACCESSING SELF CARE PRODUCTS
 Colour: **WHITE**
 Overall: **300mm (W) x 90mm (D) 235mm (H)**
 Product Weight: **0.3kg**
 Max. Load: **130kg**



Code: **RG55L**
 Desc: **SPACE SAVER SHORT LEG**
 SET OF 4 INNER LEG TUBES TO REDUCE HEIGHT OF PRODUCT BY 100mm
 Colour: **GREY**
 Overall: **Ø22.4 x 1.2mm tube**



Code: **RGSB**
 Desc: **SPACE SAVER BOWL LID**
 EASY TO INSTALL AND REMOVE
 Colour: **WHITE**
 Overall: **310mm (L) x 280mm (W) 140mm (H)**
 Product Weight: **250g BOWL 70g LID**
 Max. Load: **130kg**



Code: **RGSBLH**
 Desc: **SPACE SAVER BOWL LID HOLDER**
 POSITIONED ON EITHER SIDE OF TOILET SEAT RAISER
 Colour: **WHITE**
 Overall: **200mm (W) x 220mm (H) 50mm (D)**
 Product Weight: **0.4kg**



Code: **RGSTRH**
 Desc: **SPACE SAVER FREESTANDING TOILET ROLL HOLDER**
 IDEAL FOR EVERY COMMUNE
 Colour: **WHITE**
 Height: **810 - 860mm (2 SETTINGS)**
 Product Weight: **1.8kg**

Ⓟ PATENTED DESIGN



Code: **RG554H2**
Desc: **SPACE SAVER SHOWER STOOL**
ADJUSTABLE HEIGHT,
FLAT PACK
Colours: **WHITE**
Seat Height: **485 - 610mm**
Seat Size: **470mm (W)**
355mm (D)
Product Weight: **2.6kg**
Max. Load: **130kg**



SPACE SAVER SHOWER STOOL WITH CADDY

- Caddy clips onto handle on a:
Space Saver Shower Chair (RG555H)
Space Saver Toilet Seat Raiser (RG515)
Space Saver Shower Stool (RG554H2)
- Ideal for assisting and easily accessing self care products and toiletries such as shampoo and washing aids.
- Easily attachable with no tools required.

DEMENTIA ASSIST HANDLES SHOWN



Code: **RG5401**
Desc: **ALUM. SHOWER STOOL**
ALUM. FRAME,
FLAT PACK
ADJUSTABLE HEIGHT
Colours: **WHITE**
Seat Height: **465 - 575mm**
Seat Width: **435mm**
Product Weight: **4kg**
Max. Load: **159kg**



Code: **RG9012**
Desc: **ROUND BATH STOOL**
ALUMINIUM LEGS
ADJUSTABLE HEIGHT
Colours: **WHITE**
Seat Height: **480 - 670mm**
Seat Diameter: **Ø325mm**
Product Weight: **1.7kg**
Max. Load: **150kg**



Code: **RG9012R**
Desc: **ROUND BATH STOOL**
ROTATING SEAT,
ALUM. LEGS
ADJUSTABLE HEIGHT
Colours: **WHITE / BLUE SEAT**
Seat Height: **480 - 670mm**
Seat Diameter: **Ø325mm**
Product Weight: **1.8kg**
Max. Load: **150kg**



Ⓟ PATENTED DESIGN



Code: **RG555H**
Desc: **SPACE SAVER SHOWER CHAIR**
ADJUSTABLE HEIGHT,
FLAT PACK
Colours: **WHITE**
Seat Height: **485 - 610mm**
Seat Size: **470mm (W) x 355mm (D)**
Product Weight: **3.1kg**
Max. Load: **130kg**



Code: **RG5402**
Desc: **ALUM. SHOWER CHAIR**
ALUM. FRAME, FLAT PACK
ADJUSTABLE HEIGHT
Colours: **WHITE**
Seat Height: **445 - 545mm**
Seat Width: **460mm**
Product Weight: **5kg**
Max. Load: **159kg**



Code: **RG9400**
Desc: **FOLDING SHOWER CHAIR**
ALUMINIUM FRAME
ADJUSTABLE HEIGHT
Colours: **WHITE**
Seat Height: **460 - 530mm**
Seat Width: **470mm**
Product Weight: **3kg**
Max. Load: **100kg**



Code: **RG9095**
Desc: **DELUXE TRANSFER BENCH**
ALUMINIUM FRAME
ADJUSTABLE HEIGHT
Colours: **WHITE**
Seat Height: **450 - 570mm**
Seat Width: **800mm**
Product Weight: **4kg**
Max. Load: **130kg**



Code: **RG9201**
Desc: **BATH BOARD**
ADJUSTABLE HANDLES
PLASTIC CONSTRUCTION
Colours: **WHITE**
Seat Width: **450 - 710mm**
Seat Depth: **390mm**
Product Weight: **3kg**
Max. Load: **150kg**



Code: **RG9202**
Desc: **BATH BOARD**
ADJUSTABLE HANDLES
PLASTIC CONSTRUCTION
Colours: **WHITE**
Seat Width: **400 - 630mm**
Seat Depth: **350mm**
Product Weight: **3kg**
Max. Load: **150kg**





Code: **RG8923**
Desc: **COMMODE SHOWER TRANSPORT CHAIR**
INCLUDES SPLASHGUARD + BOWL
Colours: **WHITE**
Seat Height: **530mm**
Seat Width: **450mm**
Product Weight: **10.5kg**
Max. Load: **130kg**



Code: **RG570V**
Desc: **PLUS MOBILE SHOWER COMMODE**
INCLUDES SPLASHGUARD + BOWL
Colours: **WHITE with PALE BLUE**
Seat Height: **530mm**
Seat Width: **515mm**
Product Weight: **10.5kg**
Max. Load: **130kg**



Code: **RG570WK**
Desc: **PLUS MOBILE SHOWER COMMODE SELF PROPEL REAR WHEELS**
INCLUDES LEFT AND RIGHT HAND BRAKE
Colours: **GREY**
Overall: **Ø610mm / 24"**
Product Weight: **2.6kg (per wheel)**
Max. Load: **150kg**



Code: **RG57GP**
Desc: **PLUS MOBILE SHOWER COMMODE GEL PADS**
10mm GEL PAD TO RELIEVE PRESSURE WHEN SITTING
Colours: **BLUE**
Overall: **370mm (W) x 145mm (D) 20mm (H)**
Product Weight: **1.5kg**
Max. Load: **130kg**

NEW!

P PATENTED DESIGN



Code: **RG515**
 Desc: **SPACE SAVER TOILET SEAT RAISER**
 ADJUSTABLE HEIGHT, FLAT PACK, SPLASHGUARD
 Colours: **WHITE with SILVER LEGS**
 Seat Height: **500 - 625mm**
 Seat Width: **470mm**
 Product Weight: **3.5kg**
 Max. Load: **130kg**



Code: **RG8560**
 Desc: **FOLDING TOILET SEAT RAISER**
 ALUM. FRAME, ADJ. HEIGHT INCLUDES SPLASHGUARD + BOWL
 Colours: **WHITE**
 Seat Height: **460 - 620mm**
 Seat Width: **480mm**
 Product Weight: **3.5kg**
 Max. Load: **100kg**



Code: **RG8590**
 Desc: **TOILET SEAT RAISER**
 ALUM. FRAME, ADJ. HEIGHT INCLUDES SPLASHGUARD + BOWL
 Colours: **WHITE**
 Seat Height: **450 - 600mm**
 Seat Width: **450mm**
 Product Weight: **3.5kg**
 Max. Load: **150kg**



NEW!


R REGISTERED DESIGN




Code: **RG518**
 Desc: **FOLDING TOILET SEAT RAISER**
 ADJ. HEIGHT, INCLUDES SPLASHGUARD + BOWL
 Colours: **WHITE**
 Seat Height: **460 - 560mm**
 Seat Width: **490mm**
 Product Weight: **5kg**
 Max. Load: **130kg**

NEW!


R REGISTERED DESIGN




Code: **RG52FTS**
 Desc: **DUAL GRAB FOLDING TOILET SURROUND**
 CAN BE FOLDED FOR STORAGE OR TRANSIT
 Colours: **WHITE**
 Overall Height: **600 - 700mm**
 Internal Width: **495mm**
 Product Weight: **4kg**
 Max. Load: **130kg**



Code: **RG510-1**
 Desc: **RAISED TOILET SEAT (DONUT)**
 INCLUDES LID, PLASTIC
 Colours: **WHITE**
 Seat Height: **75mm**
 Seat Width: **360mm**
 Product Weight: **1.6kg**
 Max. Load: **100kg**




Code: **RG510-2**
 Desc: **RAISED TOILET SEAT (DONUT)**
 INCLUDES LID, PLASTIC
 Colours: **WHITE**
 Seat Height: **100mm**
 Seat Width: **360mm**
 Product Weight: **1.75kg**
 Max. Load: **100kg**



Code: **RG510-3**
 Desc: **RAISED TOILET SEAT (DONUT)**
 INCLUDES LID, PLASTIC
 Colours: **WHITE**
 Seat Height: **135mm**
 Seat Width: **360mm**
 Product Weight: **1.9kg**
 Max. Load: **100kg**

FOLDS FOR STORAGE AND TRANSIT



NEW! © REGISTERED DESIGN




Code: **RG5930**
Desc: **HANDI GRAB RAIL - 300mm**
IDEAL USE FOR TOILET, SHOWER, BATHROOM, NEAR STEPS
Colours: **WHITE**
Rail Length: **340mm**
Rail Height: **75mm**
Product Weight: **0.2kg**
Max. Load: **100kg**


NEW! © REGISTERED DESIGN



Code: **RG5945**
Desc: **HANDI GRAB RAIL - 450mm**
IDEAL USE FOR TOILET, SHOWER, BATHROOM, NEAR STEPS
Colours: **WHITE**
Rail Length: **490mm**
Rail Height: **75mm**
Product Weight: **0.3kg**
Max. Load: **100kg**



Code: **RG5905**
Desc: **MULTI ANGLE GRAB RAIL**
2 SECTION, VARIOUS ANGLES
Colours: **WHITE**
Rail Lengths: **305 / 460mm**
Rail Height: **75mm** (WALL - INSIDE OF RAIL)
Product Weight: **0.35kg / 0.45kg**
Max. Load: **227kg**



Code: **RG5906**
Desc: **MULTI ANGLE GRAB RAIL EXTENSION**
610mm EXTENSION
Colours: **WHITE**
Rail Length: **610mm**
Rail Height: **75mm** (WALL - INSIDE OF RAIL)
Product Weight: **0.65kg**
Max. Load: **227kg**



Code: **RG5903**
Desc: **SUCTION CUP GRAB RAIL - 300mm**
NO FIXINGS REQUIRED
Colours: **WHITE**
Hand Grip Length: **120mm**
Overall Length: **410mm**
Product Weight: **0.4kg**
Max. Assisted Load: **65kg**



Code: **RG5904**
Desc: **SUCTION CUP GRAB RAIL - 400mm**
NO FIXINGS REQUIRED
Colours: **WHITE**
Hand Grip Length: **220mm**
Overall Length: **510mm**
Product Weight: **0.6kg**
Max. Assisted Load: **65kg**



Code: **RG5902**
Desc: **SUCTION CUP SHOWER HEAD HOLDER**
FITS MOST SHOWER HEAD DEVICES
Colours: **WHITE**
Collar Size: **Ø25mm**
Overall Width: **60mm**
Product Weight: **70g**
Max. Load: **5kg**



MULTI ANGLE GRAB RAIL

Ideal for use in
Shower
Toilet
Bathroom
Steps

® REGISTERED DESIGN



® REGISTERED DESIGN




NEW! ® REGISTERED DESIGN



KATIEplus+
Day Chair Range

NEW! ® REGISTERED DESIGN



KATIEplus+
Day Chair Range

Code: **RG635BL / RG635GR
RG635I / RG635SG**

Desc: **LOW BACK KATIE CHAIR**

ADJUSTABLE HEIGHT,
CONTOURED PADDED
UPHOLSTERY

Colours: **BLUE / GREY / IVORY
SAGE GREEN UPHOLSTERY
GREY FRAME**

Seat Width: **480mm** (distance between armrests)

Product Weight: **9kg**

Max. Load: **200kg**

Code: **RG630BL / RG630GR
RG630I / RG630SG**

Desc: **HIGH BACK KATIE CHAIR**

ADJUSTABLE HEIGHT,
CONTOURED PADDED
UPHOLSTERY

Colours: **BLUE / GREY / IVORY
SAGE GREEN UPHOLSTERY
GREY FRAME**

Seat Width: **480mm** (distance between armrests)

Product Weight: **13kg**

Max. Load: **200kg**

Code: **RG636BL / RG636GR**
Desc: **LOW BACK KATIE PLUS
DAY CHAIR**

ADJUSTABLE HEIGHT,
CONTOURED PADDED
UPHOLSTERY

Colours: **BLUE / GREY UPHOLSTERY
GREY FRAME**

Seat Width: **530mm** (distance between armrests)

Product Weight: **13kg**

Max. Load: **200kg**

Code: **RG631BL / RG631GR**
Desc: **HIGH BACK KATIE PLUS
DAY CHAIR**

ADJUSTABLE HEIGHT,
CONTOURED PADDED
UPHOLSTERY

Colours: **BLUE / GREY UPHOLSTERY
GREY FRAME**

Seat Width: **530mm** (distance between armrests)

Product Weight: **15kg**

Max. Load: **200kg**



NEW!



Code: **RG63LRBL / RG63LRGR
RG63LRI / RG63LRSG**
Desc: **KATIE LEG REST / STOOL**
ADJUSTABLE HEIGHT+ANGLE
PADDED UPHOLSTERY
Colour: **BLUE / GREY / IVORY /
SAGE GREEN UPHOLSTERY
GREY FRAME**
Platform: **450 x 380mm**
Product Weight: **5kg**
Max. Load: **120kg** (evenly distributed)

Code: **RG63SFT**
Desc: **KATIE SIDE FOLDING TABLE**
CAN BE FITTED ON LEFT
RIGHT FOLDS DOWN WHEN
NOT REQUIRED
Colour: **BLACK**
Table Size: **400 x 200mm**
Product Weight: **3kg**
Max. Load: **10kg** (evenly distributed)

Code: **RG63RLT**
Desc: **KATIE RUBBER LEG TIPS
(SET OF 2)**
ATTACH TIPS FOR IMPROVED
GRIP ON HARD FLOORING
Colour: **GREY**
Overall: **Ø40mm (D) x 40mm (H)**
Product Weight: **0.3kg each**
Max. Load: **200kg**



**PERFECT
DAY CHAIR**

**MORE
KATIE ACCESSORIES**
(COMING SOON)
- **OXYGEN BOTTLE CARRIER**
- **I.V. POLE**

NEW!



Code: **LC0101 / LC0102 / LC0103**
 Desc: **DENVER DUAL MOTOR LIFT CHAIR**
FULL RECLINE TO STANDING
 Colour: **MOCHA / GREY / BLACK**
 Overall Width: **905mm**
 Overall Depth: **940mm**
 Overall Height: **1150mm**
 Seat Height: **490mm**
 Seat Width: **470mm**
 Product Weight: **65kg**
 Max. Load: **160kg**



DENVER LIFT CHAIR NEW FEATURES



LIFT / REST / RECLINE ARMCHAIR
 REMOVABLE BACKREST CUSHIONS
 REMOVABLE PADDING IN BACKREST CUSHIONS
 USB PORT (TYPE A) ON HANDSET
 REMOVABLE ARMREST / HEADREST COVER

LEGREST EXTENSION



Code: **LC0901**
Desc: **ULTRACARE MOBILE HEIGHT ADJUSTABLE LIFT CHAIR**
FULL RECLINE TO STANDING
Colour: **GREY**
Overall Width: **760mm**
Overall Depth: **850mm**
Overall Height: **1140mm**
Seat Height: **530 - 710mm**
Seat Width: **540mm**
Product Weight: **65kg**
Max. Load: **160kg**

ULTRACARE LIFT CHAIR

HEIGHT ADJUSTABLE / MOBILE
REMOVABLE ARMRESTS
WINGS FOR SIDE TRANSFERS
REMOVABLE LITHIUM-ION BATTERY
WITH AUDIO AND VISUAL WARNING WHEN
BATTERY REQUIRES RECHARGING



REMOVABLE
LITHIUM-ION
BATTERY





Code: **LC0301 / LC0302 / LC0303**
 Desc: **FLORENCE WALL SAVER LIFT CHAIR**
 FULL RECLINE TO STANDING
 Colours: **TOFFEE / GREY / BLACK**
 Seat Height: **480mm**
 Seat Width: **480mm**
 Product Weight: **52kg**
 Max. Load: **160kg**



Code: **LC0401 / LC0402**
 Desc: **TARANTO FABRIC 4 MOTOR LIFT CHAIR**
 FULL RECLINE TO STANDING
 Colours: **DARK GREY / TAUPE**
 Seat Height: **500mm**
 Seat Width: **470mm**
 Product Weight: **65kg**
 Max. Load: **160kg**



Code: **LCCOVER / LCCOVERGR LCCOVEBK**
 Desc: **LOUNGE CHAIR PROTECTIVE COVER**
 PROTECTS YOUR CHAIR FROM ACCIDENTAL STAINS
 Colours: **BROWN / GREY / BLACK**
 Overall: **1880mm (L) x 1590mm (W)**
 Product Weight: **0.7kg**



Code: **LCCSGY**
 Desc: **CHAIR PROTECTIVE COVER**
 PROTECTS YOUR CHAIR FROM ACCIDENTAL STAINS
 Colours: **GREY**
 Overall: **580mm (L) x 490mm (W)**
 Product Weight: **0.1kg**



Code: **RG622**
 Desc: **OVER BED CHAIR TABLE**
 TABLE TOP, FIXED BASE, RETAINING RIM
 Colour: **WHITE**
 Platform Height: **800 - 1100mm**
 Platform Size: **850mm (W) 425mm (D)**
 Product Weight: **11kg**
 Max. Load: **12kg**

Code: **RG624**
 Desc: **OVER BED CHAIR TABLE**
 ONE PIECE TOP
 Colour: **BEECH**
 Platform Height: **725 - 1130mm**
 Platform Size: **555mm (W) 325mm (D)**
 Product Weight: **12kg**
 Max. Load: **20kg**

Code: **RG625**
 Desc: **OVER BED CHAIR TABLE**
 HEIGHT ADJUSTABLE SPLIT PLATFORM
 Colour: **WHITE**
 Platform Height: **700 - 980mm**
 Platform Size: **750mm (W) 400mm (D)**
 Product Weight: **8kg**
 Max. Load: **10kg**

Code: **RG626**
 Desc: **OVER BED CHAIR TABLE**
 HEIGHT ADJUSTABLE SPLIT PLATFORM
 Colour: **BEECH**
 Platform Height: **700 - 980mm**
 Platform Size: **750mm (W) 400mm (D)**
 Product Weight: **9kg**
 Max. Load: **10kg**

Code: **RGOBCT**
 Desc: **OVER BED CHAIR TABLE**
 NEW TOP DESIGN! TILTING MAIN PLATFORM
 Colour: **BEECH**
 Platform Height: **725 - 1130mm**
 Platform Size: **555mm (W) 325mm (D)**
 Product Weight: **14kg**
 Max. Load: **20kg**



NEW!

Code: **RG4044**
 Desc: **VARI AIR FLOTATION CUSHION**
 8 x 9 CELLS, INCLUDES PUMP + COVER
 Colours: **BLACK RUBBER**
 Seat Height: **70mm**
 Cushion Size: **400mm (W) 440mm (D)**
 Product Weight: **1.6kg**
 Max. Load: **120kg**

NEW!

Code: **RG4545**
 Desc: **VARI AIR FLOTATION CUSHION**
 8 x 8 CELLS, INCLUDES PUMP + COVER
 Colours: **BLACK RUBBER**
 Seat Height: **70mm**
 Cushion Size: **450mm (W) 450mm (D)**
 Product Weight: **2.3kg**
 Max. Load: **120kg**


NEW!

Code: **RG4848**
 Desc: **VARI AIR FLOTATION CUSHION**
 8 x 8 CELLS, INCLUDES PUMP + COVER
 Colours: **BLACK RUBBER**
 Seat Height: **70mm**
 Cushion Size: **480mm (W) 480mm (D)**
 Product Weight: **2.5kg**
 Max. Load: **120kg**


Code: **RG4040B**
 Desc: **AIR CUSHION**
 8 x 8 CELLS, INCLUDES PUMP + COVER
 Colours: **BLACK RUBBER**
 Seat Height: **30mm**
 Cushion Size: **400mm (W) 420mm (D)**
 Product Weight: **0.6kg**
 Max. Load: **120kg**

PRESSURE REDISTRIBUTION

Pressure mapping proves that the **VARI AIR** Flotation Cushion redistributes pressure when sitting to provide a more comfortable seating solution



Supports the contours of the user when sitting




Code: **WC75**
 Desc: **SEAT WEDGE CUSHION**
 AUSTRALIAN MADE, TRADITIONAL FOAM, REMOVABLE WASHABLE COVER
 Colours: **BLACK**
 Cushion Height: **75mm**
 Cushion Size: **400mm (W) x 400mm (D)**
 Product Weight: **0.3kg**




Code: **WC100**
 Desc: **SEAT WEDGE CUSHION**
 AUSTRALIAN MADE, TRADITIONAL FOAM, REMOVABLE WASHABLE COVER
 Colours: **BLACK**
 Cushion Height: **100mm**
 Cushion Size: **400mm (W) x 400mm (D)**
 Product Weight: **0.4kg**




Code: **WC125**
 Desc: **SEAT WEDGE CUSHION**
 AUSTRALIAN MADE, TRADITIONAL FOAM, REMOVABLE WASHABLE COVER
 Colours: **BLACK**
 Cushion Height: **125mm**
 Cushion Size: **400mm (W) x 400mm (D)**
 Product Weight: **0.5kg**







Code: **RGCC1S**
 Desc: **PREMIUM COMFORT CUSHION SOFT**
 AUSTRALIAN MADE
 TRADITIONAL FOAM
 REMOVABLE WASHABLE COVER
 Colours: **BLACK**
 Cushion Height: **80mm**
 Cushion Size: **460mm (W) x 430mm (D)**
 Product Weight: **0.6kg**




Code: **RGCC2M**
 Desc: **PREMIUM COMFORT CUSHION MEDIUM**
 AUSTRALIAN MADE
 TRADITIONAL FOAM
 REMOVABLE WASHABLE COVER
 Colours: **BLACK**
 Cushion Height: **80mm**
 Cushion Size: **460mm (W) x 430mm (D)**
 Product Weight: **0.6kg**




Code: **RGCC3F**
 Desc: **PREMIUM COMFORT CUSHION MEDIUM / FIRM**
 AUSTRALIAN MADE
 TRADITIONAL FOAM
 REMOVABLE WASHABLE COVER
 Colours: **BLACK**
 Cushion Height: **80mm**
 Cushion Size: **460mm (W) x 430mm (D)**
 Product Weight: **0.7kg**




Code: **RGCC4T**
 Desc: **PREMIUM COMFORT CUSHION TWIN LAYER SOFT / MEDIUM**
 AUSTRALIAN MADE
 TRADITIONAL FOAM
 REMOVABLE WASHABLE COVER
 Colours: **BLACK**
 Cushion Height: **80mm**
 Cushion Size: **460mm (W) x 430mm (D)**
 Product Weight: **0.65kg**




Code: **RGCC5D**
 Desc: **PREMIUM COMFORT CUSHION DELUXE SOFT**
 AUSTRALIAN MADE
 TRADITIONAL FOAM
 REMOVABLE WASHABLE COVER
 Colours: **BLACK**
 Cushion Height: **80mm**
 Cushion Size: **460mm (W) x 430mm (D)**
 Product Weight: **0.7kg**




Code: **RGC50**
 Desc: **DUAL LAYER MEMORY FOAM COMFORT CUSHION**
 TRADITIONAL DUAL LAYER FOAM
 REMOVABLE WASHABLE COVER
 Colours: **DARK GREY**
 Cushion Height: **50mm**
 Cushion Size: **435mm (W) x 435mm (D)**
 Product Weight: **0.7kg**



Code: **RGC75**
 Desc: **DUAL LAYER MEMORY FOAM COMFORT CUSHION**
 TRADITIONAL DUAL LAYER FOAM
 REMOVABLE WASHABLE COVER
 Colours: **DARK GREY**
 Cushion Height: **75mm**
 Cushion Size: **435mm (W) x 435mm (D)**
 Product Weight: **0.8kg**



Code: **RGC100**
 Desc: **DUAL LAYER MEMORY FOAM COMFORT CUSHION**
 TRADITIONAL DUAL LAYER FOAM
 REMOVABLE WASHABLE COVER
 Colours: **DARK GREY**
 Cushion Height: **100mm**
 Cushion Size: **435mm (W) x 435mm (D)**
 Product Weight: **1.2kg**



Code: **RGCCAR**
 Desc: **SWIVEL CUSHION**
 SWIVELS 360° IN EITHER DIRECTION
 Colours: **BLUE**
 Cushion Size: **Ø430mm**
 Product Weight: **0.3kg**



Code: **RGFEPI**
 Desc: **PEDAL EXERCISER**
 DIGITAL DISPLAY FOLDS
 Colours: **RED / BLUE**
 Overall Height: **310mm**
 Overall Size: **500mm (L) x 390mm (D)**
 Product Weight: **2.4kg**



MEDICAL BED



COMPLIANT TO RELEVANT AUST / NZ STANDARDS

Code: **RG66BHL**

Desc: **COMFORT LUX HIGH LOW KING SINGLE BED**

4 SECTION MATTRESS PLATFORM, HIGH - LOW PLATFORM HEIGHT, MASSAGE

Colours: **GREY/BLACK**

Overall: **2050mm (L) x 1120mm (W) x 370mm (H)** (at highest setting)

Product Weight: **114kg**

Maximum User Weight: **200kg**

Code: **RG66BHLH**

Desc: **COMFORT LUX MEDICAL HIGH LOW KING SINGLE BED**

4 SECTION MATTRESS PLATFORM, MATTRESS RETAINER, CPR FUNCTION

Colours: **GREY/BLACK**

Overall: **2050mm (L) x 1120mm (W) x 370mm (H)** (at highest setting)

Product Weight: **116kg**

Maximum User Weight: **200kg**

Code: **RG66BS**

Desc: **COMFORT LUX STANDARD KING SINGLE BED**

5 SECTION MATTRESS PLATFORM, ANTISNORE

Colours: **GREY/BLACK**

Overall: **2030m (L) x 1050mm (W) x 350mm (H)** (at highest setting)

Product Weight: **78kg**

Maximum User Weight: **200kg**



Code: **RG66HB**

Desc: **COMFORT LUX BED HEAD**

MATCHING UPHOLSTERY TO BED SURROUND, EASY TO FIT AND REMOVE

Colours: **GREY**

Overall: **1100mm (W) x 600mm (H)**

Product Weight: **8kg**



Code: **RG66R18**

Desc: **COMFORT LUX BED RAIL 180mm**

ASSISTIVE GRAB RAIL TO AID USER IN AND OUT OF BED, ATTACH TO BED TO INCREASE STABILITY

Colours: **BLACK**

Overall: **35mm (D) x 180mm (W) x 700mm (H)**

Product Weight: **2.4kg**

Code: **RG66R35**


Desc: **COMFORT LUX BED RAIL 350mm**

ASSISTIVE GRAB RAIL TO AID USER IN AND OUT OF BED, ATTACH TO BED TO INCREASE STABILITY


Colours: **BLACK**

Overall: **35mm (D) x 350mm (W) x 700mm (H)**


Product Weight: **2.6kg**



Code: **RG66M150**
 Desc: **COMFORT LUX KING SINGLE BED MATTRESS 150mm**
 MULTI LAYER MEMORY FOAM CONSTRUCTION
 Colours: **BLUE/WHITE WITH RED STRIP**
 Overall: **2030mm (L) x 1050mm (W) x 150mm (H)**
 Product Weight: **11kg**
 Maximum User Weight: **200kg**




Code: **RG66M200**
 Desc: **COMFORT LUX KING SINGLE BED MATTRESS 200mm**
 MULTI LAYER MEMORY FOAM CONSTRUCTION
 Colours: **BLUE/WHITE WITH BLACK STRIP**
 Overall: **2030mm (L) x 1050mm (W) x 200mm (H)**
 Product Weight: **15kg**
 Maximum User Weight: **200kg**



Code: **RG66M250**
 Desc: **COMFORT LUX KING SINGLE BED MATTRESS 250mm**
 MULTI LAYER MEMORY FOAM CONSTRUCTION
 Colours: **BLUE/WHITE WITH BLUE STRIP**
 Overall: **2030mm (L) x 1050mm (W) x 250mm (H)**
 Product Weight: **19kg**
 Maximum User Weight: **200kg**



Code: **RG66WMC**
 Desc: **COMFORT LUX KING SINGLE MATTRESS INCONTINENCE COVER**
 DESIGNED TO FIT ONTO KING SINGLE MATTRESS, PROTECTS MATTRESS FROM ACCIDENTAL STAINS
 Colours: **BLUE**
 Overall: **2030mm (L) x 1070mm (W) x up to 250mm (H)**
 Product Weight: **1.6kg**



Code: **RG66CLCD**
 Desc: **COMFORT LUX COOL DREAM MEMORY FOAM PILLOW**
 PERFORATED MEMORY FOAM PILLOW BODY ALLOWS FOR INCREASED AIR FLOW, REMOVABLE OUTER COVER
 Colours: **WHITE**
 Overall: **650mm (L) x 410mm (W) x 170mm (H)**
 Product Weight: **1.3kg**



Code: **RG66CLPP**
 Desc: **COMFORT LUX PILLOW PERFECT ADJUSTABLE MEMORY FOAM PILLOW**
 ADJUSTABLE PILLOW TO SUIT INDIVIDUAL REQUIREMENTS, REMOVABLE PILLOW FILLING
 Colours: **WHITE**
 Overall: **650mm (L) x 410mm (W)**
 Overall Height: **Height and weight depends on amount of filling being used**
 Product Weight with all filling: **1.3kg**

ADJUSTABLE
 pillow filling to suit your individual requirements





Code: **RGBMT**
Desc: **BED MOVING TROLLEY**
ASSIST WITH MOVING BED / MATTRESS,
REDUCE MANUAL HANDLING
Colours: **RED**
Overall Width: **680mm**
Overall Base Height: **280mm**
Overall Height Range: **770mm - 1200mm** (platform to top strap bracket)
Product Weight: **18kg**
Maximum User Weight: **130kg**



Code: **RG64BS**
Desc: **BLANKET SUPPORT**
SUPPORTS BEDDING ABOVE PERSON IN
BED, FOLDS FOR STORAGE AND TRANSIT
Colours: **SILVER**
Overall Length: **510mm**
Overall Width: **460mm**
Overall Height: **560mm**
Product Weight: **2kg**
Max User Weight: **10kg**



Code: **RG64ABS**
Desc: **ADJUSTABLE BLANKET SUPPORT**
SUPPORTS BEDDING ABOVE PERSON IN
BED, FOLDS FOR STORAGE AND TRANSIT
Colours: **SILVER**
Overall Length: **500mm**
Overall Width: **470mm** (large frame end)
Overall Height: **565 - 690mm**
Product Weight: **2.5kg**
Max User Weight: **10kg**



Code: **RG66MC**
Desc: **MATTRESS PROTECTIVE COVER**
IDEAL FOR USE ON ANY SIZE MATTRESS,
PROTECTS MATTRESS FROM ACCIDENTAL
STAINS, CLEAN AND HYGENIC
Colours: **GREY**
Overall Length: **1980mm**
Overall Width: **895mm**
Overall Height: **1mm**
Product Weight: **0.6kg**



Code: **RG614**
Desc: **PEAK BED RAIL**
STEEP SLOPING RAIL FLAT PACK
STRONG / LIGHTWEIGHT
Colours: **SILVER**
Rail Diameter: **Ø25mm TUBE**
Rail Dims: **600mm (L) x 550mm (W) x 495mm (H)**
Product Weight: **2kg**




Code: **RG612**
Desc: **SLOPE BED RAIL**
SLOPING RAIL FLAT PACK
STRONG / LIGHTWEIGHT
Colours: **SILVER**
Rail Diameter: **Ø25mm TUBE**
Rail Dims: **600mm (L) x 550mm (W) x 495mm (H)**
Product Weight: **2kg**




Code: **RG611**
Desc: **STANDARD BED RAIL**
STRAIGHT RAIL FLAT PACK
STRONG / LIGHTWEIGHT
Colours: **SILVER**
Rail Diameter: **Ø25mm TUBE**
Rail Dims: **600mm (L) x 550mm (W) x 495mm (H)**
Product Weight: **2kg**




Code: **RG61SP**
Desc: **STORAGE POCKET**
TO SUIT RG611: STANDARD BED RAIL
IDEAL MAGAZINE / BOOK STORAGE
MESH WITH ELASTIC ENDS
Colours: **BLACK**
Pocket Dims: **550mm (L) x 200mm (H)**




Max. Height: **155mm**
Max. Angle: **15°**




Max. Height: **233mm**
Max. Angle: **15°**




Max. Height: **310mm**
Max. Angle: **15°**



Max. Height: **388mm**
Max. Angle: **15°**



Max. Height: **465mm**
Max. Angle: **15°**



Max. Height: **614mm**
Max. Angle: **17°**

Code: **RAM2**
Desc: **SINGLE FOLDING ACCESS RAMP**
ALUMINIUM, CARRY HANDLE
Colours: **SILVER**
Ramp Length: **610mm**
Ramp Width: **700mm (int)**
Product Weight: **4.3kg**
Max. Load: **270kg**


Code: **RAM3**
Desc: **SINGLE FOLDING ACCESS RAMP**
ALUMINIUM, CARRY HANDLE
Colours: **SILVER**
Ramp Length: **915mm**
Ramp Width: **700mm (int)**
Product Weight: **6.8kg**
Max. Load: **270kg**

Code: **RAM4**
Desc: **SINGLE FOLDING ACCESS RAMP**
ALUMINIUM, CARRY HANDLE
Colours: **SILVER**
Ramp Length: **1220mm**
Ramp Width: **700mm (int)**
Product Weight: **8.8kg**
Max. Load: **270kg**


Code: **RAM5**
Desc: **SINGLE FOLDING ACCESS RAMP**
ALUMINIUM, CARRY HANDLE
Colours: **SILVER**
Ramp Length: **1520mm**
Ramp Width: **700mm (int)**
Product Weight: **11.6kg**
Max. Load: **270kg**

Code: **RAM6**
Desc: **SINGLE FOLDING ACCESS RAMP**
ALUMINIUM, CARRY HANDLE
Colours: **SILVER**
Ramp Length: **1826mm**
Ramp Width: **700mm (int)**
Product Weight: **14kg**
Max. Load: **270kg**


Code: **RAM7**
Desc: **SINGLE FOLDING ACCESS RAMP**
ALUMINIUM, CARRY HANDLE
Colours: **SILVER**
Ramp Length: **2125mm**
Ramp Width: **700mm (int)**
Product Weight: **18.4kg**
Max. Load: **270kg**




Code: **RAM32**
Desc: **ONE PIECE RAMP**
ALUMINIUM, REINFORCED EDGING, NON SLIP SURFACE
Colours: **SILVER**
Ramp Length: **610mm**
Ramp Width: **705mm (int)**
Product Weight: **4.8kg**
Max. Load: **270kg**




Code: **RAM33**
Desc: **ONE PIECE RAMP**
ALUMINIUM, REINFORCED EDGING, NON SLIP SURFACE
Colours: **SILVER**
Ramp Length: **914mm**
Ramp Width: **705mm (int)**
Product Weight: **6.9kg**
Max. Load: **270kg**



Code: **RAM34**
Desc: **ONE PIECE RAMP**
ALUMINIUM, REINFORCED EDGING, NON SLIP SURFACE
Colours: **SILVER**
Ramp Length: **1217mm**
Ramp Width: **705mm (int)**
Product Weight: **8.9kg**
Max. Load: **270kg**












Code: **RAM40**
Desc: **THRESHOLD ACCESS RAMP**
ALUMINIUM NON SLIP SURFACE
Colours: **SILVER / BLACK GRIP**
Ramp Length: **330mm**
Ramp Width: **810mm**
Product Weight: **2kg**
Max. Load: **270kg**



Code: **RAM20B**
Desc: **RAMP PAIR**
ALUMINIUM NON SLIP SURFACE
Colours: **SILVER**
Ramp Length: **1500mm**
Ramp Width: **190mm**
Product Weight: **3kg each**
Max. Load: **270kg**





Code: **RAM70**
 Desc: **TELESCOPIC RAMP (PAIR)**
 ALUMINIUM, NON SLIP SURFACE
 Colour: **SILVER / BLACK GRIP**
 Ramp Length: **1060 - 1800mm**
 Ramp Width: **155mm (int)**
 Product Weight: **3.5kg each**
 Max. Load: **270kg**


Code: **RAM74T**
 Desc: **TELESCOPIC RAMP (PAIR)**
 ALUMINIUM, NON SLIP SURFACE
 Colour: **SILVER**
 Ramp Length: **740 - 1215mm**
 Ramp Width: **162mm (int)**
 Product Weight: **3.3kg each**
 Max. Load: **270kg**

Code: **RAM76T**
 Desc: **TELESCOPIC RAMP (PAIR)**
 ALUMINIUM, NON SLIP SURFACE
 Colour: **SILVER**
 Ramp Length: **1045 - 1820mm**
 Ramp Width: **162mm (int)**
 Product Weight: **4.5kg each**
 Max. Load: **270kg**

Code: **RAM78T**
 Desc: **TELESCOPIC RAMP (PAIR)**
 ALUMINIUM, NON SLIP SURFACE
 Colour: **SILVER**
 Ramp Length: **1348 - 2430mm**
 Ramp Width: **162mm (int)**
 Product Weight: **5.8kg each**
 Max. Load: **270kg**


Code: **RAM80T**
 Desc: **TRIPLE TELESCOPIC RAMP (PAIR)**
 ALUMINIUM, NON SLIP SURFACE
 Colour: **SILVER**
 Ramp Length: **1070 - 2390mm**
 Ramp Width: **162mm (int)**
 Product Weight: **6.9kg each**
 Max. Load: **270kg**

NEW!




Code: **RAMR2025**
Desc: **RUBBER WEDGE RAMP**
NON SLIP SURFACE ANGLED SIDES
Colours: **BLACK**
Wedge Height: **25mm**
Wedge Dims: **1105mm (D) x 210mm (W)**
Product Weight: **3.5kg**
Max. Load: **270kg**


NEW!




Code: **RAMR2050**
Desc: **RUBBER WEDGE RAMP**
NON SLIP SURFACE ANGLED SIDES
Colours: **BLACK**
Wedge Height: **50mm**
Wedge Dims: **1105mm (D) x 410mm (W)**
Product Weight: **19kg**
Max. Load: **270kg**



Code: **RAMR2075**
Desc: **RUBBER WEDGE RAMP**
NON SLIP SURFACE ANGLED SIDES
Colours: **BLACK**
Wedge Height: **75mm**
Wedge Dims: **1105mm (D) x 460mm (W)**
Product Weight: **20kg**
Max. Load: **270kg**




Code: **RAMR20100**
Desc: **RUBBER WEDGE RAMP**
NON SLIP SURFACE ANGLED SIDES
Colours: **BLACK**
Wedge Height: **100mm**
Wedge Dims: **1105mm (D) x 540mm (W)**
Product Weight: **21kg**
Max. Load: **270kg**




15mm

Code: **RAMR1415**
Desc: **RUBBER WEDGE RAMP**
NON SLIP SURFACE CAN BE CUT TO REQUIRED LENGTH
Colours: **BLACK**
Wedge Height: **15mm**
Wedge Dims: **900mm (D) 140mm (W)**
Product Weight: **1.6kg**
Max. Load: **270kg**




20mm

Code: **RAMR1520**
Desc: **RUBBER WEDGE RAMP**
NON SLIP SURFACE CAN BE CUT TO REQUIRED LENGTH
Colours: **BLACK**
Wedge Height: **20mm**
Wedge Dims: **900mm (D) 150mm (W)**
Product Weight: **2.1kg**
Max. Load: **270kg**




25mm

Code: **RAMR1525**
Desc: **RUBBER WEDGE RAMP**
NON SLIP SURFACE CAN BE CUT TO REQUIRED LENGTH
Colours: **BLACK**
Wedge Height: **25mm**
Wedge Dims: **900mm (D) 150mm (W)**
Product Weight: **2.5kg**
Max. Load: **270kg**




30mm

Code: **RAMR2030**
Desc: **RUBBER WEDGE RAMP**
NON SLIP SURFACE CAN BE CUT TO REQUIRED LENGTH
Colours: **BLACK**
Wedge Height: **30mm**
Wedge Dims: **900mm (L) 200mm (W)**
Product Weight: **4.2kg**
Max. Load: **270kg**




35mm

Code: **RAMR2035**
Desc: **RUBBER WEDGE RAMP**
NON SLIP SURFACE CAN BE CUT TO REQUIRED LENGTH
Colours: **BLACK**
Wedge Height: **35mm**
Wedge Dims: **900mm (D) 200mm (W)**
Product Weight: **4.7kg**
Max. Load: **270kg**



40mm

Code: **RAMR2040**
Desc: **RUBBER WEDGE RAMP**
NON SLIP SURFACE CAN BE CUT TO REQUIRED LENGTH
Colours: **BLACK**
Wedge Height: **40mm**
Wedge Dims: **900mm (D) 205mm (W)**
Product Weight: **4.9kg**
Max. Load: **270kg**



45mm

Code: **RAMR2045**
Desc: **RUBBER WEDGE RAMP**
NON SLIP SURFACE CAN BE CUT TO REQUIRED LENGTH
Colours: **BLACK**
Wedge Height: **45mm**
Wedge Dims: **900mm (D) 255mm (W)**
Product Weight: **7.2kg**
Max. Load: **270kg**

REDGUM Brand:

Registered NDIS Provider
TGA (Therapeutic Goods Administration) Registered
Proudly 100% Australian Owned and Operated
ATSA Member
Perth / Sydney Offices + Warehousing
Distributors Australia / New Zealand wide
In-house R+D team - Product Development
In-house Servicing of all products

FOR-DE Group brands include:

REDGUM Brand: Various Assisted Living and Lifestyle products

SHOPRIDER Mobility: Electric Mobility Scooters and Powerchairs

IKON Scooters: Electric Mobility Scooters


Scooter Xtras: Mobility Scooter and Powerchair Accessories

Baricare: Various Bariatric Assisted Living and Lifestyle products

REDGUM Brand owns the copyright to this product catalogue. Any information contained in this catalogue shall not be photocopied, reproduced or translated into any other languages without the prior written consent of REDGUM Brand.

REDGUM Brand strives to provide the highest quality products possible that are subject to change without notice.

© Copyright **REDGUM BRAND** 2024



REDGUM brand
ASSISTED LIVING PRODUCTS



PROUDLY
AUSTRALIAN
OWNED / OPERATED

ASK FOR THE BRAND

redgumbrand.com.au