


NEW PRODUCT



REPLACEMENT
TYRE FOR
WHEELCHAIR



CAN FIT ONTO ANY
WHEELCHAIR WITH
A TYRE SIZE

Ø24" x 1-3/8"
DIAMETER



REDGUM brand strives to provide the highest quality products possible that are subject to change without notice.

IMPORTANT NOTE:

It is the responsibility of the user to ensure the suitability of this product for its intended use. The user should consult their appropriate healthcare professional if unsure of the correct product usage for their unique circumstances.

WHEELCHAIR REPLACEMENT TYRE

650069

FEATURES

- Replacement tyre for wheelchair
- Smooth ride **FLAT FREE** tyre

Replacement tyre can fit onto the following Self-Propelled Wheelchairs:

- Titan Heavy Duty Wheelchair
(**RG22HD**)
- Onyx Tilt-in-Space Self Propel Wheelchair 18 & 20"
(**RG350SP18 & RG350SP20**)
- Pearl Deluxe Aluminium Wheelchair
(**RG18AAD & RG20AAD**)
(**RG18AADRB & RG20AARB**)
- Quartz Aluminium Wheelchair
(**RG310A18BK - BLACK**)
(**RG310A18BL - BLUE**)
(**RG310A18R - RED**)
- Crystal 2 Self-Propel Wheelchair 18 & 20"
(**RG18B & RG20B**)

SPECIFICATIONS

(All dimensions shown are nominal)

- Overall: **Ø24" x 1-3/8"**
- Product Weight: **0.6kg**

Please contact **REDGUM Brand** if you require further information on this product

P: **08 9248 4180**

E: **orders@for-de.com.au**

1 Business Way Malaga Western Australia 6090
Part of the **For-de Group Pty Ltd**

ASK FOR THE BRAND

redgumbrand.com.au