

GREAT NEW COLOUR



REDGUM brand strives to provide the highest quality products possible that are subject to change without notice.

IMPORTANT NOTE:

It is the responsibility of the user to ensure the suitability of this product for its intended use. The user should consult their appropriate healthcare professional if unsure of the correct product usage for their unique circumstances.

COMFORTLITE TRANSIT WHEELCHAIR

RG312GR18 - GREEN
RG312BL18 - STEEL BLUE

FEATURES

- Lightweight: 9.5kg (Aluminium frame)
- Fixed armrest with PU padding
- Foldable backrest with storage pocket
- Removable seat with backrest cushions
- Double Handbrake
- Fixed leg rests with flip up footrests
- Large flat free tyres
- Seat belt
- Calf strap / Side straps (both removable)
- Anti tippers

SPECIFICATIONS

(All dimensions shown are nominal)

- Overall:
 - 840mm (L) - Footrests up**
 - 610mm (W) - 270mm (Folded)**
 - 975mm (H) - 730mm (Folded)**
- Seat:
 - 480mm (H) x 460mm (W) x 420mm (D)**
- Backrest: **450mm (H) x 455mm (W)**
- Armrest: **200mm (H)**
- Wheels: **150mm / 6" (Front)**
300mm / 12" (Rear)
- Product Weight: **9.5kg**
- Max. User Weight: **100kg**

Please contact **REDGUM Brand** if you require further information on this product



P: 08 9248 4180

E: orders@for-de.com.au

1 Business Way Malaga Western Australia 6090
Part of the **For-de Group Pty Ltd**