

# HEIGHT ADJUSTABLE SEAT WALKER



REMOVABLE  
STORAGE BAG

THREE COLOURS!

RG4215

## FEATURES

- Great walking and stability aid
- Height adjustable legs and handles
- Padded backrest
- Ideal for indoor and outdoor use
- Folds up for transport and storage
- Lightweight (6.4kg)
- Included storage bag under seat
- Cushioned seat
- Dual handbrakes with park brake
- Large flat free wheels

**FOLDS UP  
FOR TRANSPORT  
& STORAGE**



TURN OVER FOR  
SPECIFICATIONS

**OCEAN BLUE**  
RG4215BL

**BURNT ORANGE**  
RG4215OG

**RUBY RED**  
RG4215RD



ASK FOR THE BRAND

[redgumbrandmobility.com.au](http://redgumbrandmobility.com.au)

# HEIGHT ADJUSTABLE SEAT WALKER



HEIGHT ADJUSTABLE  
HANDLES

HEIGHT ADJUSTABLE  
LEGS / SEAT

## SPECIFICATIONS

(All dimensions shown are nominal and are in mm)

- Overall Dimensions (at lowest setting):  
**715mm (L) x 610mm (W) x 710mm (H)**
- Folded: **715mm (L) x 435mm (W) x 710mm (H)**
- Seat (H from floor):  
**340mm (L) x 440mm (W) x 520-560mm (H)**
- Backrest (from seat to top of backrest)  
**455mm (W) x 230mm (H)**
- Handle Height (from floor to top of handle)  
Setting 1: **710mm**  
Setting 2: **735mm**  
Setting 3: **760mm**  
Setting 4: **785mm**  
Setting 5: **810mm**  
Setting 6: **835mm**  
Setting 7: **860mm**
- Wheels: **Ø200mm / 8" flat free tyres**
- Product Weight: **6.4kg**
- Max. User Weight: **120kg**

Please contact **REDGUM Brand** if you require further information on this product.

**P: 08 9248 4180**

**E: [orders@for-de.com.au](mailto:orders@for-de.com.au)**

ASK FOR THE BRAND

[redgumbrandmobility.com.au](http://redgumbrandmobility.com.au)