



IVORY
(RG635I)



BLUE
(RG635BL)

ADJUSTABLE
seat height and depth



GREY
(RG635GR)



SAGE GREEN
(RG635SG)

Backrest Height
480mm
(seat to top of backrest)

FEATURES

- Adjustable Seat Height
- Adjustable Seat Depth
- Contoured Padded Backrest
- Comfortable Padded Seat
- Large Padded Armrests
- Chair Mobility Wheels to rear legs as standard
- Backrest can be folded down for transit or storage
- Strong lightweight construction
- Epoxy Powder Coated framework
- Vinyl covered upholstery

SPECIFICATIONS

- Overall Dimensions (at lowest setting):
665mm (W) x 650mm (D) x 915mm (H)
- Seat Width: **480mm** (between armrests)
- Seat Depth: **435 - 475mm**
- Seat Height: **450 - 575mm** (from floor to top of seat)
25mm adjustment increments
- Armrest Height: **220mm** (from top of seat to top of armrest)
- Backrest Width: **480mm** (at widest point)
- Backrest Height: **480mm**
- Backrest Angle: **98°**
- Mobility Wheel: **Ø50mm**
- Product Weight: **11kg**
- Maximum User Weight: **200kg**

Please contact **REDGUM Brand** if you require further information on this product

P: 08 9248 4180
E: orders@for-de.com.au

ASK FOR THE BRAND

redgumbrandmobility.com.au