

These wheelchairs have been crash tested according to standard "ISO 7176-19" and may be used as a seat in motor vehicles.



**AVAILABLE IN
18" and 20"
SEAT WIDTH!**



REDGUM brand strives to provide the highest quality products possible that are subject to change without notice.

IMPORTANT NOTE:

It is the responsibility of the user to ensure the suitability of this product for its intended use. The user should consult their appropriate healthcare professional if unsure of the correct product usage for their unique circumstances.

ONYX - TILT IN SPACE TRANSIT WHEELCHAIRS

RG350TR18
RG350TR20

FEATURES

- Seat and backrest are angle adjustable
- Detachable armrests
- Height adjustable armrests
- Elevating, fully adjustable legrests
- Removable height adjustable headrest
- Locking brake with drum hand brake
- Adjustable anti tip caps
- Quick release rear wheels
- Solid tyres
- Pressure reduction padded seat, legs, arms, backrest and headrest support

SPECIFICATIONS

- Overall length: **1050 - 1090mm** (with footrests)
Overall length: **830 - 870mm** (without footrests)
Overall width: **580 - 730mm**
Overall width: **650 - 700mm** (without armrest)
Overall height **990 - 1060mm** (without headrest)
Overall height: **710 - 780mm**
(with folded backrest, without armrests)
Seat width: **450mm and 500mm**
Seat depth: **460 - 520mm**
Seat height: **460 - 530mm** (with cushion)
Backrest height: **550 - 585mm**
Backrest angle: **0° - 54°**
Front Wheels: **8" - 200mm**
Rear Wheels: **24" - 600mm**
Rear Wheels: **16" - 400mm**
Distance between footrest and seat: **400 - 570mm**
Foot plate angle: **80° - 100°**
Angle between seat and footrest: **104° - 194°**
Static stability downhill: **10°**
Static stability uphill: **15°**
Static stability sideways: **12°**
Obstacle climbing: **60mm**
Product Weight: **38kg**
Maximum Weight Loading: **135kg**

Please contact **REDGUM Brand** if you require further information on this product

P: 08 9248 4180
E: orders@for-de.com.au

1 Business Way Malaga Western Australia 6090
Part of the **For-de Group Pty Ltd**

ASK FOR THE BRAND

redgumbrand.com.au