

**NEW PRODUCT!**

**Convenient  
Carrier**



**Fits on  
Most Mobility  
Solutions**



# PADDED OXYGEN BOTTLE CARRIER

SE312

## FEATURES

- Supportive padding to protect your oxygen bottle
- Attaches to most mobility solutions
- Adjustable straps
- Easy to use buckles
- Included drawstring keeps carrier closed

## SPECIFICATIONS


(All specifications may vary +/- 10mm)

- Suits **B, C and D** medical gas cylinders
- **Overall:** 590mm (H) x  $\varnothing$ 200mm
- **Product Weight:** 0.4kg
- **Maximum Load:** 20kg

**Drawstring  
Keeps Carrier  
Closed**



**Adjustable  
Straps with Easy  
to Use Buckle**



**P: 08 9248 4180**  
**E: [orders@for-de.com.au](mailto:orders@for-de.com.au)**



For more information  
about this product

**ASK FOR THE BRAND**

**scooterxtras.com.au**