



IMPROVE YOUR ACCESS !



Raised non slip surface

REDGUM brand strives to provide the highest quality products possible that are subject to change without notice.



**Improve access for
mobility equipment**

IMPORTANT NOTE:

It is the responsibility of the user to ensure the suitability of this product for its intended use.

The user should consult their appropriate healthcare professional if unsure of the correct product usage for their unique circumstances.

RUBBER WEDGE RAMPS

FEATURES

- Solid black rubber construction
- Raised non slip surface provides excellent grip
- Suitable for indoor and outdoor use
- Can be cut down or grooved to suit application
- Maximum weight capacity: **270kg**
- Available in various sizes - see below:

RAMR1415

900mm (W) x 140mm (D) x 15mm (H)
Weight: **1.6kg**

RAMR1520

900mm (W) x 150mm (D) x 20mm (H)
Weight: **2.1kg**

RAMR1525

900mm (W) x 150mm (D) x 25mm (H)
Weight: **2.5kg**

RAMR2030

900mm (W) x 200mm (D) x 30mm (H)
Weight: **4.2kg**

RAMR2035

900mm (W) x 200mm (D) x 35mm (H)
Weight: **4.7kg**

RAMR2040

900mm (W) x 200mm (D) x 40mm (H)
Weight: **4.9kg**

RAMR2045

900mm (W) x 250mm (D) x 45mm (H)
Weight: **7.2kg**

Please contact **REDGUM** if you require any more information on this product.

P: 08 9248 4180

E: orders@for-de.com.au

REDGUM brand - 1 Business Way, Malaga, Western Australia 6090
Part of the **For-de Group Pty Ltd**

ASK FOR THE BRAND

redgumbrand.com.au