

NEW SEAT WIDTH  
**490mm**



**QUICK RELEASE  
SELF PROPEL  
WHEELS**

**SWING AWAY  
LEG RESTS**

# CRYSTAL<sup>2</sup> 20" SELF PROPEL WHEELCHAIR

RG20B

## FEATURES

- Lightweight Aluminium frame
- Compact folding frame and fold down backrest
- Removable / height adjustable swing away leg rests
- Fold up footplates with heel straps
- Attendant operated hand / park brakes
- Occupant operated park brakes
- Padded armrests
- Quick release self propel rear wheels

## SPECIFICATIONS

(All dimensions shown are nominal)

- Overall: **730mm W x 1080mm L x 925mm H**  
(including swing away leg rests)
- Overall Folded: **350mm W x 670mm L x 735mm H**  
(Leg rests / rear wheels removed)
- Seat: **490mm W x 420mm D x 520mm H**
- Armrest: **300mm L x 50mm W x 220mm H**
- Backrest: **550mm W x 400mm H**
- Footrest: **175mm W x 155mm D x 340 - 440mm H**
- Wheels: **Ø200mm / 8" (Front) - Ø610mm / 24" (Rear)**
- Main Frame Weight: **8.9kg**
- Wheel Weight: **1.9kg each**
- Leg Rest Weight: **0.6kg each**
- Max. User Weight: **127kg**

Please contact **REDGUM Brand** if you require further information on this product

**COMPACT  
FOLDING  
FRAME**



P: **08 9248 4180**  
E: **orders@for-de.com.au**

1 Business Way Malaga Western Australia 6090  
Part of the **For-de Group Pty Ltd**

**IMPORTANT NOTE:**

It is the responsibility of the user to ensure the suitability of this product for its intended use. The user should consult their appropriate healthcare professional if unsure of the correct product usage for their unique circumstances.