

ANTI SLIP TIPS!



MADE OF STRONG DURABLE RUBBER

PROVIDES STABILITY AND GRIP

RG16
ø16mm



RG19
ø19mm



RG26
ø24.8mm



RG22
ø22mm



RG23
ø25.4mm

RG27
ø27.4mm

RG28
ø28.6mm

ANTI SLIP & ANTI SKID RUBBER BASE

IMPORTANT NOTE:

It is the responsibility of the user to ensure the suitability of this product for its intended use. The user should consult their appropriate healthcare professional if unsure of the correct product usage for their unique circumstances.

WALKING STICK TIPS

RG16 | RG19 | RG22 | RG23 | RG26
RG27 | RG28

FEATURES

- Inner metal ring - lengthens tip life
- Provides stability and grip
- Anti slip and anti skid rubber base
- Made of strong and durable rubber
- Suction grip assures weight placement and easy release

SPECIFICATIONS

RG16

To suit: ø16mm shaft
Footprint: ø30mm
Max. Load: 130kg
Colours: Black/Grey

RG19

To suit: ø19mm shaft
Footprint: ø37.5mm
Max. Load: 130kg
Colours: Black

RG22

To suit: ø22mm shaft
Footprint: ø45mm
Max. Load: 130kg
Colours: Grey

RG23

To suit: ø25.4mm shaft
Footprint: ø40mm
Max. Load: 130kg
Colours: Grey

RG26

To suit: ø24.8mm shaft
Footprint: ø57.5mm
Max. Load: 130kg
Colours: Grey

RG27

To suit: ø27.4mm shaft
Footprint: ø40mm
Max. Load: 130kg
Colours: Grey

RG28

To suit: ø28.6mm shaft
Footprint: ø40mm
Max. Load: 130kg
Colours: Grey

Please contact REDGUM Brand if you require further information on this product.

P: 08 9248 4180

E: orders@for-de.com.au

REDGUM brand: 1 Business way Malaga Western Australia 6090
Part of the For-de Group Pty Ltd

redgumbrand.com.au