

HEIGHT ADJUSTABLE !



ADJUSTABLE HEIGHT SEAT WALKER

RG4210ADJ - BLUE | RG4210ADJ - RED

FEATURES

- Height adjustable seat and handles
- User friendly loop brakes
- Ergonomic handgrips to reduce pressure
- Comfortable padded seat
- Padded backrest loop
- Storage bag
- Large diameter wheels for excellent mobility
- Flat free rubber tyres
- Swivelling front wheels
- Strong lightweight aluminium frame
- Folds up for transit or storage

SPECIFICATIONS

Overall Dimensions (At lowest setting):

600mm (W) x 760mm (H)

660mm (Depth - Unfolded), 300mm (Folded)

Handle height (From floor to top of handle):

Low: 760 - 900mm (7 settings)

High: 830 - 970mm (7 settings)

Seat:

440mm (W - Between handles) x 330mm (D)

540 / 590mm (From floor to top of seat)

Wheels:

8" / Ø200mm (Front - Swivel / Rear - Fixed)

Product weight:

7.5kg

Max. user weight:

120kg

Please contact
REDGUM if you
require any more
information on
this product.



P: 08 9248 4180

E: orders@for-de.com.au

REDGUM brand - 1 Business way, Malaga, Western Australia 6090
Part of the FOR-DE Group Pty Ltd.

redgumbrand.com.au

IMPORTANT NOTE:

It is the responsibility of the user to ensure the suitability of this product for its intended use.

The user should consult their appropriate healthcare professional if unsure of the correct product usage for their unique circumstances.