

REDGUM brand

ASSISTED LIVING PRODUCTS

NEW COLOUR!

**SUPER
LIGHTWEIGHT**
(Aluminium frame)

**EASY TO
ASSEMBLE**
(no tools required)



RG4402:
TRI WALKER
CHAMPAGNE



RG4402BL:
TRI WALKER
SKY BLUE



RG4402PR:
TRI WALKER
PURPLE



TRI WALKER

CHAMPAGNE (RG4402)
SKY BLUE (RG4402BL)
PURPLE (RG4402PR)

FEATURES

- **SUPER LIGHTWEIGHT!**
- Epoxy powder coated **aluminium** frame
- Modern contemporary design
- Large wheels for improved mobility
- Removable zip lock storage bag with carry handle
- Height adjustable handles
- Locking handbrakes
- Easy to engage / release brake levers
- Ergonomic rubber hand grips
- Swivel front wheel
- Easy to assemble - no tools required
- Folds up for storage and transport

SPECIFICATIONS

(All dimensions shown are nominal)

- Overall: **580mm (L) x 675mm (W) x 830mm (H)**
(lowest setting)
- Handle Heights: **830 - 930mm: 5 settings**
(floor to top of handle)
- Handle Width: **480mm** (overall)
- Wheels: **190mm Diameter**
- Product Weight: **4.6kg** (without bag)
- Max User Weight: **115kg**

Please contact **REDGUM Brand** if you require further information on this product.

P: **08 9248 4180**
E: **orders@for-de.com.au**

ASK FOR THE BRAND

redgumbrandmobility.com.au

IMPORTANT NOTE: It is the responsibility of the user to ensure the suitability of this product for its intended use. The user should consult their appropriate healthcare professional if unsure of the correct product usage for their unique circumstances.

REDGUM brand strives to provide the highest quality products possible that are subject to change without notice.