




SKY BLUE
RG4415BL

PURPLE
RG4415PR

Two
Colour
Options



Stands
when
folded



LOW SEAT SIDE FOLDING WALKER

FEATURES

- Great walking and stability aid
- Height adjustable handles
- Ergonomic hand grips to reduce pressure
- Brake and park brakes
- Convenient built in seat
- Lightweight aluminium construction
- Walking stick holder fitted as standard
- Removable storage bag fitted as standard
- Large diameter wheels for excellent mobility
- Flat free tyres
- Walker can stand up right when folded
- Folds up easily for transit or storage

SPECIFICATIONS (All dims shown are nominal)

- Overall Width: **615mm (260mm folded)**
- Overall Length: **730mm**
- Handle Height: **885 - 985mm (5 settings)**
- Seat: **460mm (W) x 250mm (D)**
- Backrest Depth: **360mm**
- Seat Height: **510mm (floor to top of seat)**
- Wheel Diameter: **Ø200mm / 8"**
- Product Weight: **7.5kg**
- Max. User Weight: **136kg**

Please contact **REDGUM Brand** if you require further information on this product

P: **08 9248 4180**
E: **orders@for-de.com.au**

ASK FOR THE BRAND

redgumbrandmobility.com.au