

NEW PRODUCT



EXCELLENT  
WALKING AID

**IMPORTANT NOTE:**

It is the responsibility of the user to ensure the suitability of this product for its intended use.

The user should consult their appropriate healthcare professional if unsure of the correct product usage for their unique circumstances.

REDGUM brand strives to provide the highest quality products possible that are subject to change without notice.

# HEAVY DUTY QUAD CANE WITH OFFSET HANDLE

RG5611

## FEATURES

- Excellent walking / balancing aid
- Suitable for left and right handed users
- Height adjustable with locking collar
- Lightweight (Aluminium / Steel construction)
- Large low profile offset quad base
- Soft foam hand grip
- Rubber tips

## SPECIFICATIONS

Top Handle Height: **740-960mm**  
(25mm increment settings)


Base: L: **305mm** W: **210mm** H: **80mm**

Max User Weight: **227kg**

Product Weight: **1.6kg**

Please contact **REDGUM Brand**  
if you require further information on  
this product

HEIGHT ADJUSTABLE  
WITH SOFT FOAM  
HAND GRIPS



P: **08 9248 4180**

E: **orders@for-de.com.au**

1 Business Way Malaga Western Australia 6090  
Part of the **For-de Group Pty Ltd**

ASK FOR THE BRAND

redgumbrand.com.au