

**HEIGHT  
ADJUSTABLE!**



### FEATURES

- Plastic platform with retaining lip
- Easy to clean, sealed platform surface
- Positive locking adjustable height platform
- Large height adjustment lever
- Practical "U" shaped base frame to fit around chairs
- Low clearance base to fit under beds
- Twin wheel locking castors
- Beige epoxy powder coated metal frame
- Chrome mild steel upright

### SPECIFICATIONS

Platform Dimensions:

**W: 850mm x D: 425mm x H: 800 - 1100mm**

Base Frame:

**W: 640mm x H: 100m**

Castors:  $\phi$ 50mm - plastic twin wheel locking castors

Product weight: 11kg

Max. Load: 12kg - evenly distributed

Please contact **REDGUM Brand** if you require further information on this product



#### IMPORTANT NOTE:

It is the responsibility of the user to ensure the suitability of this product for its intended use. The user should consult their appropriate healthcare professional if unsure of the correct product usage for their unique circumstances.

**P: 08 9248 4180**

**E: [orders@for-de.com.au](mailto:orders@for-de.com.au)**

1 Business way Malaga Western Australia 6090  
Part of the **For-de Group Pty Ltd**