

TILT FUNCTION ON PLATFORM!



TWIN WHEEL
LOCKING CASTORS



OVER BED/ CHAIR TABLE

RG623

FEATURES

- Vinyl wrapped platform with retaining lip
- Height adjustable platform
- Easy to adjust tilt function on platform
- Easily collapses for transport or storage
- Low clearance base to fit under beds
- Twin wheel locking castors
- Grey epoxy powder coated metal frame
- Chrome mild steel upright

SPECIFICATIONS

All dimensions shown are nominal and at lowest setting

Platform:

600mm (W) x 500mm (D) x 685 - 965mm (H)

Platform angles: **0° / 15° / 30° / 45° / 90°**

Base Frame:

Internal Width: **610mm**

Frame Height: (Under bed clearance): **100mm**

Flattened:


1020mm (L) x 630mm (W) x 120mm (H)

Castors: **ø50mm**

Product Weight: **9.5kg**

Maximum Load: **12kg**

Please contact
REDGUM Brand
if you require
further information
on this product.



IMPORTANT NOTE:

It is the responsibility of the user to ensure the suitability of this product for its intended use. The user should consult their appropriate healthcare professional if unsure of the correct product usage for their unique circumstances.

P: 08 9248 4180

E: orders@for-de.com.au

REDGUM brand: 1 Business way Malaga Western Australia 6090
Part of the For-de Group Pty Ltd

redgumbrand.com.au