

HEIGHT ADJUSTABLE  
PLATFORM!



TWIN WHEEL  
LOCKING CASTORS



**IMPORTANT NOTE:**

It is the responsibility of the user to ensure the suitability of this product for its intended use. The user should consult their appropriate healthcare professional if unsure of the correct product usage for their unique circumstances.

# OVER BED/ CHAIR TABLE

RG624

## FEATURES

- Stable working surface
- Easy to clean surface - vinyl wrapped
- Height adjustable platform
- Low clearance base to fit under beds
- Twin wheel locking castors
- Epoxy powder coated metal frame
- Chrome mild steel upright

## SPECIFICATIONS

(All dimensions shown are nominal and at lowest setting)

Platform:

**805mm (W) x 385mm (D) x 715 - 1125mm (H)**

Base Frame:


Internal Width: **800mm**

Frame Height: (Under bed clearance): **90mm**

Castors: **ø50mm**

Product Weight: **11kg**

Maximum Load: **12kg**



HEIGHT  
ADJUSTABLE  
PLATFORM

Please contact **REDGUM Brand** if you require further information on this product.

P: **08 9248 4180**

E: **orders@for-de.com.au**

REDGUM brand: **1 Business Way Malaga Western Australia 6090**  
Part of the **For-de Group Pty Ltd**