

NEW PRODUCT



Neatly  
folds away



Height  
adjustable



Removable  
wrist strap

Ergonomic hand grip  
with narrow neck for  
improved comfort



For more information  
about this product

# CARBON FIBRE

## WALKING STICK

RGWS390

### FEATURES

- Durable carbon fibre body
- Height adjustable
- Ergonomic grip handle
- Anti rattle locking device
- Convenient removable wrist strap
- Large 19mm rubber non-slip tip
- Collapse for storage and transport

### SPECIFICATIONS

All dimensions shown are nominal and at lowest setting

Length:

**760mm - 860mm (285mm folded)**

Handle: **125mm x 23mm**

Height Settings:

**760mm, 785mm, 810mm,  
835mm, 860mm**

Product Weight: **250g**

Max. User Weight: **113kg**

P: **08 9248 4180**

E: **orders@for-de.com.au**

ASK FOR THE BRAND

redgumbrandmobility.com.au