

# CROSS COUNTRY SEAT WALKER

RG4416

**Feature**  
Packed Design



**Adjustable**  
Handle Height



**Sleek**  
Profile



**Folds for**  
Storage



**Large 12"**  
Front Wheels



Walking Stick & Bottle Holder  
**Included**



## FEATURES

- Designed to travel across uneven ground
- Large flat free wheels for excellent mobility & stability
- Thick padded contoured backrest
- Adjustable handle height
- Walking stick holder included
- Bottle holder included
- Removable carry bag
- Hand and park brakes
- Folds for storage and transport

## SPECIFICATIONS

(All specifications may vary +/- 10mm)

- **Overall (at lowest setting):**  
790mm (L) x 640mm (W) x 890mm (H)
- **Seat Dimensions:** 220mm (L) x 480mm (W) x 600mm (H)
- **Handlebar Height:** 860mm - 930mm
- **Wheel Diameter:** 300mm / 12" (Front)  
240mm / 8" (Rear)
- **Product Weight:** 10.2kg
- **Maximum User Weight:** 136kg

P: **08 9248 4180**  
E: **orders@for-de.com.au**



For more information  
about this product

ASK FOR THE BRAND

[redgumbrandmobility.com.au](http://redgumbrandmobility.com.au)